

Lymphastim

BTL THERAPEUTIC ENCYCLOPAEDIA



1	Beauty/Aesthetic/Wellness	3
1.1	Cellulite	3
1.2	Detoxification.....	4
1.3	Heavy legs syndrome.....	5
1.4	Improvement in skin trophic	6
1.5	Lipoedema	7
1.6	Postliposuction care	8
1.7	Reconvalescence	9
1.8	Regeneration.....	10
1.9	Tenalgia - lower limbs	11
1.10	Tenalgia - upper limbs.....	12
1.11	Varicose veins prevention	13
2	Medical	14
2.1	Acrocyanosis - lower limbs.....	14
2.2	Acrocyanosis - upper limbs	15
2.3	Algodystrophic syndrome - lower limbs.....	16
2.4	Algodystrophic syndrome - upper limbs	17
2.5	Arthrosis - lower limbs.....	18
2.6	Arthrosis - upper limbs	19
2.7	Carpal Tunnel Syndrome - upper limbs.....	20
2.8	Diabetic angiopathy.....	21
2.9	Gout, holarthritis.....	22
2.10	Chronic oedemas	23
2.11	Ischemic disease of lower limbs.....	24
2.12	Neuropathy, Neuralgia	25
2.13	Obesity - lower limbs	26
2.14	Obesity - upper limbs	27
2.15	Oncology	28
2.16	Posttraumatic care - lower limbs	29
2.17	Posttraumatic care - upper limbs.....	30
2.18	Raynaud's syndrome - upper limbs	31
2.19	Repair of poorly performed MLD	32
2.20	Trophic disorders of skin - lower limbs	33
2.21	Trophic disorders of skin - upper limbs.....	34
2.22	Venous insufficiency.....	35



1 BEAUTY/AESTHETIC/WELLNESS

1.1 CELLULITE

BEAUTY/AESTHETIC/WELLNESS

<i>Program</i>	Y-1019 program: Y-1019 Cellulite
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 5:00 ----- 2. program: 0002 time: 15:00 ----- 3. program: 0009 time: 40:00
<i>Recommended accessories</i>	arm, leg , trousers
<i>Length of application</i>	60 min
<i>Frequency of treatments</i>	min. twice a week
<i>Number of treatments</i>	as needed
<i>Pressure</i>	40 - 80 mmHg
<i>Effect</i>	analgesic, antiedematous, increase of lymphatic drainage, support of venous return, myorelaxation, support of microcirculation, support of tissue regeneration, support of nutrition of the skin and subcutis, trophic, increase of local metabolism (1544)



1.2 DETOXIFICATION

<i>Program</i>	Y-1023 program: Y-1023 Detoxification
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 10:00 ----- 2. program: 0004 time: 10:00 ----- 3. program: 0013 time: 40:00
<i>Recommended accessories</i>	arm, leg , trousers
<i>Length of application</i>	60 min
<i>Frequency of treatments</i>	min. twice a week
<i>Number of treatments</i>	4 - 6
<i>Pressure</i>	30 - 70 mmHg
<i>Effect</i>	analgesic, antiedematous, increase of lymphatic drainage, support of venous return, myorelaxation, support of microcirculation, support of tissue regeneration, increase of local metabolism (1548)



1.3 HEAVY LEGS SYNDROME

BEAUTY/AESTHETIC/WELLNESS

<i>Program</i>	Y-1022 program: Y-1022 Heavy legs syndrome
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 5:00 ----- 2. program: 0002 time: 10:00 ----- 3. program: 0009 time: 30:00
<i>Recommended accessories</i>	leg , trousers
<i>Length of application</i>	45 min
<i>Frequency of treatments</i>	min. once a week
<i>Number of treatments</i>	as needed
<i>Pressure</i>	40 - 70 mmHg
<i>Effect</i>	increase of lymphatic drainage, support of venous return, increase of lactate metabolism, support of microcirculation, support of tissue regeneration, myorelaxation (1547)



1.4 IMPROVEMENT IN SKIN TROPHIC

BEAUTY/AESTHETIC/WELLNESS

<i>Program</i>	Y-1024 program: Y-1024 Improvement in skin trophic
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 5:00 ----- 2. program: 0009 time: 25:00
<i>Recommended accessories</i>	arm, leg , trousers
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	min. twice a week
<i>Number of treatments</i>	as needed
<i>Pressure</i>	30 - 45 mmHg, with higher pressure the risk of skin rupture increases
<i>Effect</i>	analgesic, antiedematous, support of microcirculation, support of tissue regeneration, henotic, support of nutrition of the skin and subcutis, trophic
<i>Note</i>	May be used as effective part of the treatment in dermatology, especially in wound management. (1549)



1.5 LIPOEDEMA

BEAUTY/AESTHETIC/WELLNESS

<i>Program</i>	Y-1021 program: Y-1021 Lipoedema
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 35 mmHg ----- 1. program: 0003 time: 10:00 ----- 2. program: 0007 time: 20:00 ----- 3. program: 0015 time: 30:00
<i>Recommended accessories</i>	arm, leg , trousers
<i>Length of application</i>	60 min
<i>Frequency of treatments</i>	min. twice a week
<i>Number of treatments</i>	as needed
<i>Pressure</i>	20 - 50 mmHg
<i>Effect</i>	increase of lymphatic drainage, support of venous return, increase of fat metabolism, support of microcirculation, support of tissue regeneration, support of nutrition of the skin and subcutis
<i>Note</i>	Regular application improves the quality and appearance of the skin. (1546)



1.6 POSTLIPOSUCTION CARE

BEAUTY/AESTHETIC/WELLNESS

<i>Program</i>	Y-1025 program: Y-1025 Postliposuction care
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 5:00 ----- 2. program: 0002 time: 10:00 ----- 3. program: 0007 time: 30:00
<i>Recommended accessories</i>	arm, leg , trousers
<i>Length of application</i>	45 min
<i>Frequency of treatments</i>	min. twice a week
<i>Number of treatments</i>	as needed
<i>Pressure</i>	40 - 70 mmHg
<i>Effect</i>	analgesic, antiedematous, increase of lymphatic drainage, support of venous return, support of microcirculation, support of tissue regeneration, henotic, support of nutrition of the skin and subcutis, trophic
<i>Note</i>	Early application after liposuction may decrease edemas in operated area and prevent unevenness of the surface. (1550)



1.7 RECONVALESCENCE

BEAUTY/AESTHETIC/WELLNESS

<i>Program</i>	Y-1026 program: Y-1026 Reconvalence
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 35 mmHg ----- 1. program: 0006 time: 10:00 ----- 2. program: 0010 time: 10:00 ----- 3. program: 0005 time: 25:00
<i>Recommended accessories</i>	arm, leg , trousers
<i>Length of application</i>	45 min
<i>Frequency of treatments</i>	min. twice a week
<i>Number of treatments</i>	as needed, until repair
<i>Pressure</i>	30 - 70 mmHg
<i>Effect</i>	analgesic, antiedematous, myorelaxation, support of microcirculation, support of tissue regeneration, henotic
<i>Effect</i>	The pressure should be set up to patient's threshold, comfortable pressure. (1551)



1.8 REGENERATION

BEAUTY/AESTHETIC/WELLNESS

<i>Program</i>	Y-1029 program: Y-1029 Regeneration
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0002 time: 10:00 ----- 2. program: 0004 time: 10:00 ----- 3. program: 0009 time: 25:00
<i>Recommended accessories</i>	arm, leg , trousers
<i>Length of application</i>	45 min
<i>Frequency of treatments</i>	min. twice a week
<i>Number of treatments</i>	as necessary
<i>Pressure</i>	40 - 90 mmHg
<i>Effect</i>	increase of lymphatic drainage, support of venous return, increase of lactate metabolism, support of microcirculation, support of tissue regeneration, myorelaxation
<i>Note</i>	The pressure should be set up to patient's threshold, comfortable pressure. (1554)



1.9 TENALGIA - LOWER LIMBS

BEAUTY/AESTHETIC/WELLNESS

<i>Program</i>	Y-1027 program: Y-1027 Tenalgia - lower limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 40 mmHg ----- 1. program: 0014 time: 5:00 ----- 2. program: 0002 time: 15:00 ----- 3. program: 0005 time: 10:00
<i>Recommended accessories</i>	leg , trousers
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	min. twice a week
<i>Number of treatments</i>	min. 6 treatments
<i>Pressure</i>	30 - 70 mmHg
<i>Effect</i>	analgesic, antiedematous, myorelaxation, support of microcirculation, support of regeneration of tissues, henotic
<i>Note</i>	The pressure should be set up to patient's threshold, comfortable pressure. (1552)



1.10 TENALGIA - UPPER LIMBS

BEAUTY/AESTHETIC/WELLNESS

<i>Program</i>	Y-1028 program: Y-1028 Tenalgia - upper limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 40 mmHg ----- 1. program: 0014 time: 5:00 ----- 2. program: 0002 time: 15:00 ----- 3. program: 0005 time: 10:00
<i>Recommended accessories</i>	arm
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	min. twice a week
<i>Number of treatments</i>	min. 6 treatments
<i>Pressure</i>	30 - 70 mmHg
<i>Effect</i>	analgesic, antiedematous, myorelaxation, support of microcirculation, support of regeneration of tissues, henotic
<i>Note</i>	The pressure should be set up to patient's threshold, comfortable pressure. (1553)



1.11 VARICOSE VEINS PREVENTION

BEAUTY/AESTHETIC/WELLNESS

<i>Program</i>	Y-1020 program: Y-1020 Varicose veins prevention
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0002 time: 10:00 ----- 2. program: 0015 time: 20:00
<i>Recommended accessories</i>	leg , trousers
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	twice a week
<i>Number of treatments</i>	as needed
<i>Pressure</i>	30 - 50 mmHg
<i>Effect</i>	support of microcirculation, alleviation of pain caused by local ischaemia, support of neovascularisation, support of venous return
<i>Note</i>	Do not apply on acute inflammation of veins known or suspected!! (1545)



2 MEDICAL

2.1 ACROCYANOSIS - LOWER LIMBS

MEDICAL

<i>Program</i>	Y-1003 program: Y-1003 Acrocyanosis - lower limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 35 mmHg ----- 1. program: 0004 time: 5:00 ----- 2. program: 0005 time: 10:00 ----- 3. program: 0015 time: 15:00
<i>Recommended accessories</i>	leg , trousers
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	1-2 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	20 - 35 mmHg
<i>Effect</i>	analgesic, support of microcirculation, neovascularisation, support of venous return, support of tissue regeneration, henotic
<i>Note</i>	Higher pressure may increase the risk of capillary rupture! (1527)



2.2 ACROCYANOSIS - UPPER LIMBS

MEDICAL

<i>Program</i>	Y-1004 program: Y-1004 Acrocyanosis - upper limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 35 mmHg ----- 1. program: 0004 time: 5:00 ----- 2. program: 0005 time: 10:00 ----- 3. program: 0015 time: 15:00
<i>Recommended accessories</i>	arm
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	1-2 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	20 - 35 mmHg
<i>Effect</i>	analgesic, support of microcirculation, neovascularisation, support of venous return, support of tissue regeneration, henotic
<i>Note</i>	Higher pressure may increase the risk of capillary rupture! (1528)



2.3 ALGODISTROPHIC SYNDROME - LOWER LIMBS

MEDICAL

<i>Program</i>	Y-1005 program: Y-1005 Algodistrophic syndrome - lower limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 5:00 ----- 2. program: 0002 time: 10:00 ----- 3. program: 0015 time: 15:00
<i>Recommended accessories</i>	leg , trousers
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	1 - 2 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	30 - 60 mmHg
<i>Effect</i>	analgesic, antiedematous, myorelaxation, support of microcirculation, support of tissue regeneration, henotic (1529)



2.4 ALGODISTROPHIC SYNDROME - UPPER LIMBS

MEDICAL

<i>Program</i>	Y-1006 program: Y-1006 Algodistrophic syndrome - upper limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 5:00 ----- 2. program: 0002 time: 10:00 ----- 3. program: 0017 time: 15:00
<i>Recommended accessories</i>	arm
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	1 - 2 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	30 - 60 mmHg
<i>Effect</i>	analgesic, antiedematous, myorelaxation, support of microcirculation, support of tissue regeneration, henotic (1530)



2.5 ARTHROSIS - LOWER LIMBS

MEDICAL

<i>Program</i>	Y-1008 program: Y-1008 Arthrosis - lower limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 5:00 ----- 2. program: 0002 time: 15:00 ----- 3. program: 0011 time: 10:00
<i>Recommended accessories</i>	leg , trousers
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	1 - 2 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	30 - 100 mmHg
<i>Effect</i>	analgesic, antiedematous, myorelaxation, support of microcirculation, support of tissue regeneration, henotic
<i>Note</i>	Respect the feeling of the patient, adjust the pressure and gradient according to the threshold level! (1532)



2.6 ARTHROSIS - UPPER LIMBS

MEDICAL

<i>Program</i>	Y-1009 program: Y-1009 Arthrosis - upper limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 5:00 ----- 2. program: 0002 time: 15:00 ----- 3. program: 0017 time: 10:00
<i>Recommended accessories</i>	arm
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	1 - 2 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	30 - 100 mmHg
<i>Effect</i>	analgesic, antiedematous, myorelaxation, support of microcirculation, support of tissue regeneration, henotic
<i>Note</i>	Respect the feeling of the patient, adjust the pressure and gradient according to the threshold level! (1533)



2.7 CARPAL TUNNEL SYNDROME - UPPER LIMBS

MEDICAL

<i>Program</i>	Y-1000 program: Y-1000 Carpal Tunnel Syndrome - upper limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0017 time: 5:00 ----- 2. program: 0002 time: 25:00
<i>Recommended accessories</i>	arm
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	1 - 2 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	30- 60 mmHg
<i>Effect</i>	analgesic, antiedematous, myorelaxation, support of microcirculation, support of tissue regeneration, henotic (1522)



2.8 DIABETIC ANGIOPATHY

MEDICAL

<i>Program</i>	Y-0015 program: Y-0015 Diabetic angiopathy
<i>Therapy parameters</i>	protocol: Lymphatic drainage The number of filling up table: 15
<i>Recommended accessories</i>	arm, leg , trousers
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	1 - 2 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	20 - 35 mmHg
<i>Effect</i>	analgesic, support of microcirculation, neovascularisation, support of venous return, support of tissue regeneration, henotic
<i>Note</i>	Higher pressure may increase the risk of capillary rupture! (1526)



2.9 GOUT, HOLAARTHITIS

MEDICAL

<i>Program</i>	Y-1010 program: Y-1010 Gout, holarthritis
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 3:00 ----- 2. program: 0002 time: 20:00 ----- 3. program: 0014 time: 7:00
<i>Recommended accessories</i>	arm, leg , trousers
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	1-2 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	30 -100 mmHg
<i>Effect</i>	analgesic, antiedematous, myorelaxation, support of microcirculation, support of tissue regeneration, henotic
<i>Note</i>	Respect the feeling of the patient, adjust the pressure and gradient according to the threshold level! (1534)



2.10 CHRONIC OEDEMAS

MEDICAL

<i>Program</i>	Y-1007 program: Y-1007 Chronic oedemas
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 10:00 ----- 2. program: 0007 time: 20:00 ----- 3. program: 0013 time: 30:00
<i>Recommended accessories</i>	arm, leg , trousers
<i>Length of application</i>	min. 60 min
<i>Frequency of treatments</i>	min. twice a week
<i>Number of treatments</i>	as needed
<i>Pressure</i>	30 - 50 mmHg
<i>Effect</i>	antiedematous, regulation of lymphatic flow, support of lymphatic drainage, support of nutrition of the skin and subcutis, trophic, support of microcirculation, support of tissue regeneration, analgesic
<i>Note</i>	Do not exceed 70 mmHg pressure! Higher pressure may cause the rupture of lymphatic veins and the congestion of fluid in interstitium. Always respect the feeling and needs of the patient! Manual unblock of the main lymphatic nodes is recommended. (1531)



2.11 ISCHEMIC DISEASE OF LOWER LIMBS

MEDICAL

<i>Program</i>	Y-0015 program: Y-0015 Ischemic disease of lower limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage The number of filling up table: 15
<i>Recommended accessories</i>	leg , trousers
<i>Length of application</i>	45 min
<i>Frequency of treatments</i>	1 - 2 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	20 - 45 mmHg
<i>Effect</i>	analgesic, support of microcirculation, neovascularisation, support of venous return, support of tissue regeneration, henotic
<i>Note</i>	Higher pressure may increase the risk of capillary rupture! (1525)



2.12 NEUROPATHY, NEURALGIA

MEDICAL

<i>Program</i>	Y-1002 program: Y-1002 Neuropathy, Neuralgia
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 35 mmHg ----- 1. program: 0004 time: 15:00 ----- 2. program: 0002 time: 15:00
<i>Recommended accessories</i>	arm, leg , trousers
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	1 - 2 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	30 - 50 mmHg
<i>Effect</i>	analgesic, antiedematous, myorelaxation, support of neural regeneration, support of microcirculation, support of tissue regeneration, support of nutrition of the skin and subcutis, trophic, increase of local metabolism
<i>Note</i>	In acute neuropathy do not exceed the recommended pressure, respect the feeling of the patient! (1524)



2.13 OBESITY - LOWER LIMBS

MEDICAL

<i>Program</i>	Y-1011 program: Y-1011 Obesity - lower limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 5:00 ----- 2. program: 0003 time: 10:00 ----- 3. program: 0009 time: 15:00
<i>Recommended accessories</i>	leg , trousers
<i>Length of application</i>	min. 30 min
<i>Frequency of treatments</i>	min. twice a week
<i>Number of treatments</i>	as needed
<i>Pressure</i>	40 - 100 mmHg
<i>Effect</i>	increase of lymphatic drainage, support of venous return, increase of fat metabolism, support of microcirculation, support of tissue regeneration, myorelaxation, support of nutrition of the skin and subcutis
<i>Note</i>	Respect the feeling of the patient, adjust the pressure and gradient according to the threshold level! (1535)



2.14 OBESITY - UPPER LIMBS

MEDICAL

<i>Program</i>	Y-1012 program: Y-1012 Obesity - upper limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 5:00 ----- 2. program: 0003 time: 10:00 ----- 3. program: 0009 time: 15:00
<i>Recommended accessories</i>	arm
<i>Length of application</i>	min. 30 min
<i>Frequency of treatments</i>	min. twice a week
<i>Number of treatments</i>	as needed
<i>Pressure</i>	40 - 100 mmHg
<i>Effect</i>	increase of lymphatic drainage, support of venous return, increase of fat metabolism, support of microcirculation, support of tissue regeneration, myorelaxation, support of nutrition of the skin and subcutis
<i>Note</i>	Respect the feeling of the patient, adjust the pressure and gradient according to the threshold level! (1536)



2.15 ONCOLOGY

MEDICAL

<i>Program</i>	Y-1015 program: Y-1015 Oncology
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 5:00 ----- 2. program: 0007 time: 25:00 ----- 3. program: 0013 time: 15:00
<i>Recommended accessories</i>	arm, leg , trousers
<i>Length of application</i>	min. 45 min
<i>Frequency of treatments</i>	2 - 3 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	30 - 50 mmHg
<i>Effect</i>	analgesic, antiedematous, regulation of lymphatic flow, support of lymphatic drainage, support of nutrition of the skin and subcutis, trophic, henotic
<i>Note</i>	Do not exceed 70 mmHg pressure! Higher pressure may cause the rupture of lymphatic veins and the congestion of fluid in interstitium. Always respect the feeling and needs of the patient! Manual deblocation of the main lymphatic nodes is recommended. (1539)



2.16 POSTTRAUMATIC CARE - LOWER LIMBS

MEDICAL

<i>Program</i>	Y-1016 program: Y-1016 Posttraumatic care - lower limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0010 time: 5:00 ----- 2. program: 0005 time: 15:00 ----- 3. program: 0001 time: 10:00
<i>Recommended accessories</i>	leg , trousers
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	min. once a week
<i>Number of treatments</i>	as needed, until repair
<i>Pressure</i>	30 - 50 mmHg
<i>Effect</i>	analgesic, antiedematous, myorelaxation, increase of lymphatic drainage, support of venous return, support of microcirculation, support of tissue regeneration, henotic, support of nutrition of the skin and subcutis, trophic
<i>Note</i>	Respect the feeling of the patient. (1540)



2.17 POSTTRAUMATIC CARE - UPPER LIMBS

MEDICAL

<i>Program</i>	Y-1017 program: Y-1017 Posttraumatic care - upper limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0016 time: 5:00 ----- 2. program: 0005 time: 15:00 ----- 3. program: 0001 time: 10:00
<i>Recommended accessories</i>	arm
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	min. once a week
<i>Number of treatments</i>	as needed, until repair
<i>Pressure</i>	30 - 50 mmHg
<i>Effect</i>	analgesic, antiedematous, myorelaxation, support of microcirculation, support of tissue regeneration, henotic, support of nutrition of the skin and subcutis, trophic
<i>Note</i>	Respect the feeling of the patient. (1541)



2.18 RAYNAUD'S SYNDROME - UPPER LIMBS

MEDICAL

<i>Program</i>	Y-1001 program: Y-1001 Raynaud's syndrome - upper limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 20 mmHg ----- 1. program: 0015 time: 20:00 ----- 2. program: 0005 time: 10:00
<i>Recommended accessories</i>	arm
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	1 - 2 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	20 - 35 mmHg
<i>Effect</i>	analgesic, support of microcirculation, neovascularisation, support of venous return, support of tissue regeneration, henotic
<i>Note</i>	Higher pressure may increase the risk of capillary rupture! (1523)



2.19 REPAIR OF POORLY PERFORMED MLD

MEDICAL

<i>Program</i>	Y-1018 program: Y-1018 Repair of poorly performed MLD
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 30 mmHg ----- 1. program: 0001 time: 15:00 ----- 2. program: 0007 time: 45:00
<i>Recommended accessories</i>	arm, leg , trousers
<i>Length of application</i>	min. 60 min
<i>Frequency of treatments</i>	2 -3 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	30 - 50 mmHg
<i>Effect</i>	antiedematous, regulation of lymphatic flow, support of lymphatic drainage, support of nutrition of the skin and subcutis, trophic, improvement in vegetative reaction caused by poorly performed MLD, unblock of main lymphatic nodes
<i>Note</i>	Do not exceed 70 mmHg pressure! Higher pressure may cause the rupture of lymphatic veins and the congestion of fluid in interstitium. (1542)



2.20 TROPHIC DISORDERS OF SKIN - LOWER LIMBS

MEDICAL

<i>Program</i>	Y-1013 program: Y-1013 Trophic disorders of skin - lower limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 35 mmHg ----- 1. program: 0006 time: 5:00 ----- 2. program: 0011 time: 5:00 ----- 3. program: 0015 time: 20:00
<i>Recommended accessories</i>	leg , trousers
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	min. once a week
<i>Number of treatments</i>	as needed
<i>Pressure</i>	30 - 45 mmHg
<i>Effect</i>	analgesic, antiedematous, support of microcirculation, support of tissue regeneration, henotic, support of nutrition of the skin and subcutis, neovascularisation, trophic
<i>Note</i>	Higher pressure may cause the skin integrity problem! (1537)



2.21 TROPHIC DISORDERS OF SKIN - UPPER LIMBS

MEDICAL

<i>Program</i>	Y-1014 program: Y-1014 Trophic disorders of skin - upper limbs
<i>Therapy parameters</i>	protocol: Lymphatic drainage sequence pause: no recommended pressure: 35 mmHg ----- 1. program: 0006 time: 5:00 ----- 2. program: 0017 time: 5:00 ----- 3. program: 0015 time: 20:00
<i>Recommended accessories</i>	arm
<i>Length of application</i>	30 min
<i>Frequency of treatments</i>	min. once a week
<i>Number of treatments</i>	as needed
<i>Pressure</i>	30 - 45 mmHg
<i>Effect</i>	analgesic, antiedematous, support of microcirculation, support of tissue regeneration, henotic, support of nutrition of the skin and subcutis, neovascularisation, trophic
<i>Note</i>	Higher pressure may cause the rupture of the skin! (1538)



2.22 VENOUS INSUFFICIENCY

MEDICAL

<i>Program</i>	Y-0015 program: Y-0015 Venous insufficiency
<i>Therapy parameters</i>	protocol: Lymphatic drainage The number of filling up table: 15
<i>Recommended accessories</i>	arm, leg , trousers
<i>Length of application</i>	45 min
<i>Frequency of treatments</i>	1 - 2 weekly
<i>Number of treatments</i>	as needed
<i>Pressure</i>	20 - 45 mmHg
<i>Effect</i>	analgesic, support of microcirculation, neovascularisation, support of venous return, support of tissue regeneration, henotic
<i>Note</i>	Do not apply on acute inflammation of veins known or suspected !! (1521)

