



## C521M **RECUMBENT** CYCLE

Semi-recumbent positioning, extreme comfort and many pedaling options make the C521M a perfect choice for multiple patient types, especially for those with knee and hip ROM limitations.

- Bi-directional pedaling for better flexibility and rehab variations
- Adjustable cranks to vary pedal rotation for one or both sides
- Reclining seat back and pivoting arm rests for ultimate fit and comfort
- Cool mesh seat back with lumbar pouch for heat and cold pack therapy
- Ultra low start up resistance for infirm or deconditioned users

<b>USER WEIGHT</b>	<b>DIMENSIONS</b>	<b>OVERALL WEIGHT</b>	<b>POWER REQUIREMENT</b>	<b>RESISTANCE LEVELS</b>
500 lbs / 227 kg	67 x 26 x 50 in / 170 x 66 x 127 cm	183 lbs / 83 kg	15A, 120V AC, grounded	20 by 0.1 increments

### Readouts

Workout Level, Calories, Time, RPM, Distance, METS, WATTS, Speed, Heart Rate, Target Heart Rates

### Warranty

Full Commercial  
**Frame** Lifetime  
**Parts** 5 Years\*  
**Labor** 3 Year\*  
**Wear Items** 2 Years

### Power Requirements

15A, 120V AC, grounded

### Workout Programs

Hill, Custom Interval, Manual/Track, Random, Cardio/Weight Loss HRC, Custom HR, Iso Kinetic, Iso Strength, 4 User IDs

### Console



Tri-Color dot matrix display

### Optional



Medical Pedal



\*U.S. and Canada only.