

Professional Products

Clinical supplies and institutional rehab & fitness equipment



.. Mettler Sonicator Plus 930 Combination ..

Sonicator Plus 930

Modèle: OC-3272

Ultrasound and Stim Combination (Interferential, Premod, Russian)

Introducing a hot new combo that performs muscle stimulation and therapeutic ultrasound. Sonicator® Plus 930 Combo features a 1- and 3-MHz 5 sq. cm applicator for therapeutic ultrasound, plus the three most popular stim waveforms – Interferential, Premodulated and Medium Frequency (Russian).

In one compact package, the Sonicator® Plus 930 gives you:

- Independent control of up to three treatment protocols
- Two intensity controls and two intensity displays, enabling the operator to monitor all aspects of a treatment simultaneously
- An easy-to-use membrane panel designed for quick set-up and simple operation
- A convenient recessed cradle that holds the ultrasound applicator between treatments

Canadian Medical Device License No.1188

Specifications:

Input: 90-240V AC, 50-60 Hz

Weight: 8.7 pounds

Dimensions: 6"H x 12"W x 12"D

Warranty: 2 year limited

Ultrasound Specifications:

Effective Radiating Area: 5 sq. cm.

Frequency: 1.0 MHz, 3.2MHz

Modes: Continuous, Pulsed, 20%

Max. Output Power: 11W

Max. Intensity: 2.2 W/sq.cm.

Stimulation Specifications:

Output: Interferential, 0-65mA RMS, 1 kOhm load; Premodulated, 0-50 mA RMS, 1 kOhm load; Medium Frequency, 0-50 mA RMS, 1 kOhm load

Frequency: Interferential, 4000 Hz, 0-250 Hz

Output Options: Continuous, Vector Rotation, Surge, Reciprocation