

Product Specification Sheet



www.endorphin.net

370 Series Upper / Lower Platform



resistance system

Ships via *Freight*
Shipping Dimensions

US / Metric
320 lbs. / 146 kg
L - 60" / 152 cm
W - 31" / 79 cm
H - 65" / 165 cm

Product Dimensions
237 lbs. / 108 kg
L - 68" / 173 cm
W - 28" / 71 cm
H - 43" / 109 cm

specifications

Vendor Part Number: 370-E1
Suggested Retail Price: \$2895.00
Power Requirements: None
Options: 400 Multi-Positional Hand Grip
426 Foot Plates
Cover: New millennium cover design eliminates fly-wheel pinch points for improved user safety.
Flywheel: Bi-directional flywheel on the e₁ allows the user to exercise with resistance in both the forward and reverse directions.
Assembly: Partial assembly required.

Displays: Equipped with a multi-function electronic pulse meter, which reads distance, time, speed, calories, and pulse.
Resistance: Adjustable magnetic resistance features 8 incremental lever positions, which allows the user to have consistent resistance changes.
Handles: Comfort grip handles are provided to allow the user to rotate his arms forward and reverse.
Operation: The e-system resistance modules are operated and controlled independently and give the user flexibility in customizing treatment for specific injuries or for a general exercise program.

Adjustment: The height of the resistance system crank adjusts from 36" to 50" from the floor.
Positioning: The semi-recumbent seat adjusts front to rear for proper user positioning and rotates 360° to provide easy access. The seat assembly may be quickly removed to allow for wheelchair accessibility.



Endorphin's industry best warranty:

- 1 year bumper to bumper on all parts
- 5 years on all moving parts
- 20 years on the frame



6901 90th Ave. North, Pinellas Park, FL 33782
phone 727-545-9848 toll free 800-940-9844 fax 727-546-0613
www.endorphin.net sales@endorphin.net

