


EconoPrima™



The EconoPrima massage table is a Canadian made table, whisper quiet and built for patient comfort. It features an adjustable face rest platform and crescent cushion designed for proper positioning and relaxation. The arm rest and face rest platforms can be removed to accommodate treatment techniques that require the therapist to position themselves closely to the patient. It complies with the Canadian electrical safety standards.

1 SECTION TABLE

Electric Hi-Lo Control	Foot pedal
Height Adjustment	20.5" - 37.5"
Head Section	Adjustable crescent bracket with cushion
Arm Board	✓
Frame Color	Grey
Upholstery Foam	3" thick
Antimicrobial Naugahyde Upholstery	Black



SPECIFICATIONS

Dynamic Lifting Capacity	350 lbs.
Dimensions	29" x 76.75"
Power Supply	120 Volt AC 60 Hz. Max. 2.5 Amp.
Fire Retardant Upholstery	✓
CSA Equivalent	✓

Warranty

5 Years on Frame, 3 Years on Power Supply, Actuator, 1 Year on Upholstery and Foot Pedal