


ATHLETIC TRAINING EQUIPMENT CATALOG



PROTEAM STRONG - Over 5,000 installations and building...



Pro-FORM™ Vinyl Colors (CAL133 – Boston Fire Rated tops are available at an up charge. Specify "BF" when ordering)



ABOVE PRINTED COLORS MAY VARY FROM ACTUAL MATERIAL

Hausmann Industries is pleased to present the latest edition of the PROTEAM[™] catalog, representing the full line of athletic training room furnishings and products focused on the medical and rehabilitation markets. This catalog features innovative designs and product enhancements that meet the evolving roles and changing settings of today's athletic trainer – such as Hausmann's tablet-based facility design program using the MagicPlan[™] app (see page 7).

For over 15 years, PROTEAM has been the choice of athletic trainers for their training room furnishings. PROTEAM products feature durable, maintenance-free laminate surfaces for years of worry-free use. Our rugged Pro-Form vinyl with two-inch thick high-density foam tops ensures the extra durability needed in high-volume settings. Standard 3mm edging on all taping stations, treatment tables, carts, and cabinetry provides unparalleled durability.

Our modular taping stations fit any space and need. With exceptional versatility, modules can be added, changed, or relocated as necessary. PROTEAM offers a wide range of treatment tables, from traditional open storage tables to cabinet-style tables featuring optional backrests and elevating split legs. Add embossed or debossed logos to any taping station, treatment table or stool to personalize athletic training rooms. Logos are now available on mat platforms, featuring our new PROPRINT^M technology.

PROTEAM is the market leader - from taping stations and treatment tables to stools, carts, and cabinets - backed by over 60 years of manufacturing expertise.

PROTEAM STRONG - Over 5,000 installations and building...

David Hausmann

Chief Executive Officer Hausmann Industries, Inc.

George Batchelor

Dir of Sales & Marketing Hausmann Industries, Inc.

Joe Patten

National Sales Manager PROTEAM



We are proud to manufacture all of our wood products in the U.S.A. since 1955.



Athletic Training Equipment

TABLE OF CONTENTS

	3-6	7 MagicPlan
New Product Showcase	Modular Taping Stations	
8-9 RAGIN CAJUNS		11-12
Logos	Treatment and Taping Tables	Cabinet Tables
13-15 Split Leg Tables	17-19 Specialty Treatment Tables	20-21 S&W by Hausmann
22-25 Power Treatment and Exam Tables	26-29 Mat Platforms	30-37
39-40 Stools	41-43	45-48 Cabinets, Lockers, and Wardrobes

These colorful icons will appear throughout our catalog as a quick, easy-to-recognize reference to product features:

Industry Leading Quick Ship Program



The program continues to improve: more products, shorter lead times, addition of our new 1 Day QS and On-Time Performance of over 98%. Products are guaranteed to ship from factory in 1, 3, 5, or 10 business days after receipt of approved order or we pay the freight. Additional transit time

varies by destination and carrier. Quick Ship guarantee is for a maximum of 3 items.



"Peace of Mind" 3-Year Limited Warranty



Weight Capacity

400 The maximum weight this product can accommodate. Weight capacity is shown next to each product.

Bariatric Equipment


This product is designed to accommodate obese individuals. Product weight capacity is 500 lbs. or greater.

Wheelchair Accessible

This product is designed for wheelchair access.

UL Listed ίų

This product is listed by Underwriters Laboratory and meets UL 60601 Medical Standards for Safety. All UL Listed products meet Canadian UL (CUL).



3mm Rounded edges on door and drawer fronts will not dent or peel for enhanced safety and durability.



Green-Line Product

Eco-friendly construction for safer indoor environment.

NEW PRODUCTS



MODULAR TAPING STATIONS

COMMON FEATURES

PROTEAM MODULAR TAPING STATIONS:

 Upholstered Taping Modules are 36"H or 42"H.
 Storage Modules with laminate tops are 4" lower.



450

- You can order individual taping or storage modules. See page 6 or go to www.proteamtables.com.
- Durable dowel and glue laminated wood construction for long term structural integrity.
- Self-closing laminate drawer with nylon rollers.
- Doors have concealed, adjustable, self-closing 110° Euro-hinges.
- All doors and drawer fronts have 3mm high impact edging.
- · Cabinet backs are finished.
- Taping station end units have side open storage with (1) adjustable shelf.
- Comfortable 2" thick high-density foam upholstered tops are removable without tools.

COLORS: (SEE INSIDE FRONT COVER) Laminate* - Natural Oak, Wild Cherry, or Folkstone Gray

Upholstery* – 10 Pro-Form[™] vinyl colors *Custom laminate and upholstery available



A. MODEL A9521 2-Seat taping station with options

SHOWN L TO R:

- · A918 Full Extension Pullout
- · A926 Fixed 75° Upholstered Back
- A923 Set of (3) Drawers
- · A924 Waste Hamper
- A945 Top Drawer Lock
- A944 Full Height Doors
- A946 Lock/Bolt

A929 Side Pull-Out Stand					
MODEL LENGTH DEPTH HEIGHT					
A9521	6'	42"	36"		

Reverse side has solid back

B. MODEL A9400

Shown as left hand unit (standard configuration)

Option A926 Fixed 75° Upholstered Back

OPTIONS MODELS A9400, A9401

*A926 - Fixed 75° Upholstered Back *A945 - Door/Drawer Lock *99-R - Pullout, drawers, and door are on left when facing taping station

*FACTORY INSTALLED



B. MODEL A9400 Single Seat taping station

- Full extension pullout for tape storage opens to the side.
- Two 8" self-closing laminate drawers with nylon rollers for extra deep storage to the side for easy access to supplies.
- Large side storage compartment with shelf includes door with concealed, adjustable, self-closing 110° Euro-hinges.
- · Left hand configuration is standard.
- Taping station can be configured left or right. Please specify for right hand configuration

MODEL	WIDTH	DEPTH	HEIGHT
A9400	28"	42"	36"
A9401	28"	42"	42"

MODULAR TAPING STATIONS

A. MODEL A952X 2 Seat Modular Taping Station — 6 Ft. Long

MODEL	HEIGHT	DEPTH	
A9520 *	36"	36"	
A9521	36"	42"	
A9522	36"	48"	
A9523*	42"	36"	
A9524	42"	42"	
A9525	42"	48"	

*A926 not available on 36" depth.

B. MODEL A953X 3 Seat Modular Taping Station — 10 FT. Long

MODEL	HEIGHT	DEPTH	
A9530*	36"	36"	
A9531	36"	42"	
A9532	36"	48"	
A9533*	42"	36"	
A9534	42"	42"	
A9535	42"	48"	

*A926 not available on 36" depth.

C. MODEL A954X 4 Seat Modular Taping Station — 14 Ft Long

MODEL	HEIGHT	DEPTH		
A9540*	36"	36"		
A9541	36"	42"		
A9542	36"	48"		
A9543*	42"	36"		
A9544	42"	42"		
A9545	42"	48"		
*A926 not available on 36" depth.				

Visit www.proteamtables.com/literature.html for Taping Station Design Templates









D. MODEL A9520-346-24 Quick-Ship 2-Seat Taping Station With Hamper



SPECIAL FEATURES:

- Same as A9520 but includes A924 Hamper. No other options available.
- This Quick-Ship item is Oak laminate in size shown with choice of 10 Pro-Form[™] vinyl colors.

001015.			
MODEL	LENGTH	DEPTH	HEIGHT
A9520-346-24	6'	36"	36"

See Common Features on page 3.

MODULAR TAPING STATIONS





FULL EXTENSION PULL-OUT





A944* Full Height Doors




A923*
(3) DRAWERS (in lieu of doors)



A941* LAMINATE TOP WITH 3MM EDGING




VELCRO WALL PAD, 24"W x 24"H Shown mounted to wall



WASTE HAMPER AND CUTOUT




A946* LOCK-BOLT A947 4" SATIN FINISH WIRE PULL



HINGED SIDE ACCESS DOOR




A927* FRONT PULL-OUT STAND AND DRAWER Shown with Option A918



FIXED 75° UPHOLSTERED BACK, 26"H x 10"D Only available on 42"D or 48"D units





NOTE: All dimensions shown are nominal. Please define exact length if wall to wall.



camera or draw room dimensions add structural attributes

equipment into your room plan

- Add laminate & vinyl options and accessories to your **PROTEAM™** products
- Create your own objects to customize your plan
- Add rooms to your plan
- Update and modify your plan as often as needed
- Export your plan in PDF, JPEG, DXF & HTML formats


For more information visit: www.proteamtables.com/roomdesign

or email us: design@proteamtables.com









Athletic Training Room layouts created using the MagicPlan[™] App featuring PROTEAM[™] products

INSTALLATIONS












1.1







Note: the use of installation images and logos throughout this catalog does not imply any endorsement on the part of the respective institutions.



SHOW YOUR TEAM PRIDE... Make your training room POP!

Brand your equipment with your team logo embossed or debossed onto our heavy duty vinyl upholstery.

It's PROTEAM Tough & Cleanable - unlike silk screening, decals or embroidery.

OPTION	DESCRIPTION
A900*	6" Deboss Plate (Stool Tops only)
A901*	12" Deboss Plate
A902	Deboss Logo
A903*	12" Emboss Plate
A904	Emboss Logo ONE Color
A905	Emboss Logo TWO Color
A906	Emboss Logo THREE Color
A907	Emboss Logo FOUR Color
A908	Emboss Logo FIVE Color
-L	Debossed Stool Logo (Stool Tops only)
* Vector art	is preferred for proper logo processing

* Vector art is preferred for proper logo processing.

Contact your PROTEAM dealer for custom logo ordering information.

Scan QR code for artwork specifications.





Your team logo graphics are permanently embedded into the high pressure laminate surface... so logo will never peel or fade away. Washable, wear resistant and carefree.

Logo can be printed on vertical 24"W x 30"H back panel and/or on 24"W storage cabinet horizontal work surface.

Standard laminate background color is Black HPL, however other background color choices can be used & will match the (10) Pro-Form[™] upholstery colors shown on inside front cover.

Contact your dealer for laminate logo ordering information at time of order. Scan the QR code above for artwork specifications.

OPTION	DESCRIPTION	
A910*	Laminate Logo Setup Fee	
A911	Laminate Logo Wall Plate	
A912	Laminate Logo Storage Top	

* Vector art is preferred for proper logo processing. Your team logo graphics can be embedded on sides, door and drawer fronts (See page 16).

TREATMENT TABLES/TAPING TABLES



COMMON FEATURES

TREATMENT TABLES

MODELS A9002, A9024, A9041, A9087:

- 30" x 78" upholstered top with comfortable 2" thick high-density foam.
- · High-pressure laminate legs, apron and H-Brace.
- · Enhanced leg hardware system.
- · Legs have protective black vinyl scuff caps.

COLORS: (SEE INSIDE FRONT COVER) Laminate - Natural Oak or Wild Cherry

Upholstery – 10 Pro-Form[™] vinyl colors

A. MODEL A9087 600 LB. TRAINERS TABLE

SPECIAL FEATURES:

• (6) Legs and plywood top for added strength. Laminate shelf with (2) long stretchers. .

MODEL	LENGTH	WIDTH	HEIGHT
A9087	78″	30"	31"

NPTINNS

MODELS A9002, A9024, A9041, A9087

*69 – Gas Spring Backrest 0° to 75° 🜆

*FACTORY INSTALLED

B. MODEL A9024





B. MODELS A9002, A9024, A9041 TREATMENT TABLE

SPECIAL FEATURES:

- 400
- · A9024 has full width laminate shelf. A00/1 bas laminate shalf and enen

·	A9041 has laminate shelf and open storage
	cabinet.
·	Additional sizes available.

MODEL	LENGTH	WIDTH	HEIGHT	
A9024	78"	30"	31"	
A9002 (not shown) same as A9024 w/o shelf				
A9041 (inset) same as A9024 w/open storage				

C. MODEL A9093 H-BRACE TAPING TABLE

SPECIAL FEATURE:

 36" is convenient height for tapir 	ng.
--	-----

MODEL	LENGTH	WIDTH	HEIGHT
A9093	48"	27"	36"



C. MODEL A9093

D. MODEL A9097 H-BRACE TAPING TABLE WITH END CABINET

SPECIAL FEATURE:

 End facing open storage cabinet. 				
MODEL	LENGTH	WIDTH	HEIGHT	
A9097	48"	27"	36"	

NEW 27" WIDTH



ENHANCED LEG HARDWARE SYSTEM





metal threaded inserts for durable connection. Easily tighten and re-tighten as needed.

COMMON FEATURES TAPING TABLES -

MODELS A9093, A9097:



- H-Brace is recessed into table leg with unique 4-sided "Lock-Tite" joint.
- High pressure laminate legs, apron and stretchers.
- Legs have protective black vinyl scuff . caps.
- Comfortable 2" thick high-density foam . upholstered top.

COLORS: (SEE INSIDE FRONT COVER) Laminate - Natural Oak Upholstery – 10 Pro-Form[™] vinyl colors

A. MODEL A9066 TREATMENT TABLE WITH OPEN SHELF STORAGE



- Open shelf storage, accessible from both sides of the table.
- · Center section has fixed shelves. Left and right sections have (1) adjustable shelf.
- · Comfortable 2" thick high-density foam upholstered tops.
- Natural Oak, Wild Cherry, or Folkstone Gray laminate.
- · 10 Pro-Form[™] vinyl colors (specify).

MODEL	LENGTH	WIDTH	HEIGHT
A9066	78″	30"	31"

B. MODEL A9067 CABINET TREATMENT TABLE



- (6) Hinged doors on both sides of table.
- Storage accessible from both sides of table. . Center section has fixed shelves. Left and right •
- sections have (1) adjustable shelf.
- · Comfortable 2" thick high-density foam upholstered tops.
- Natural Oak, Wild Cherry, or Folkstone Gray laminate.

 10 Pro-Form[™] vinyl colors (specify). 				
MODEL	LENGTH	WIDTH	HEIGHT	
A9067	78″	30"	31"	







C. MODEL A9073 CABINET TABLE

*FACTORY INSTALLED

OPTIONS



to head section.)

- · (3) Two-door cabinets with laminate interiors and adjustable shelf.
- · Solid finished back.
- · Comfortable 2" thick high-density foam upholstered tops.
- Natural Oak, Wild Cherry, or Folkstone Gray laminate.

· 10 Pro-Form™ vinyl colors (specify).

MODEL	LENGTH	WIDTH	HEIGHT
A9073	78"	30"	31"





A. MODEL 4839 ORTHOPEDIC/HAND SURGERY TABLE New Lockable slide-out shelf. Shelf is 24"W x 13"D and fully adjustable in depth with locking knob.

- Two fully enclosed door cabinets.
- Each door cabinet is equipped with (1) adjustable laminate shelf.
- · Concealed, self-closing Euro-hinges.
- · Two all laminate drawers on nylon rollers.
- · 2" High-density urethane foam top.
- Natural Oak, Wild Cherry, or Folkstone Gray laminate.

Choice of (10) Pro-Form[™] vinyl colors.

MODEL	LENGTH	WIDTH	HEIGHT
4839	78"	30"	31"



OPTIONS

- MODEL 4079, 4143, 4839
- *07 Nose Cutout
- 23 Paper Dispenser
- 25 Paper Cutter
- 27 Paper Dispenser & Cutter Combo32 Pillow, full size
- **35** Pillow, small
- 55 Pillow, small
- *69 Gas Spring Backrest 0° to 75° (400 (4079 only)
- *99L Left Hand Table Model
- *A946 (1) Lock/Bolt for Double Doors

*FACTORY INSTALLED

- B. MODEL 4143 SPACE-SAVER CABINET TREATMENT TABLE
- Space-Saver Tables meet your needs for a full size table in a compact space.
- Full width 20" long foot end folds down to reduce table length to 52".
- Positive action locking device allows the foot end to be easily elevated into three positions: 55°, 75°, and 90°.
- Fully enclosed all laminate cabinet.
- Two hinged doors with 3mm edging.
- Concealed, self-closing Euro-hinges.
- · Laminate interiors with (1) adjustable shelf.
- · Black laminated plywood base.
- Natural Oak, Wild Cherry, or Folkstone Gray laminate.
- · Choice of (10) Pro-Form[™] vinyl colors.

MODEL	LENGTH	WIDTH	HEIGHT
4143	52" - 72"*	24"	31"

* Fully extended

C. MODEL 4079 Modality treatment table



- Designed to save time and effort with modality equipment at your fingertips.
- Spacious upholstered top with 2" high-density urethane foam.
- Two fully enclosed door cabinets with adjustable shelf.
- Mid-section of table has full extension 10"H x 23 ¾"W equipment compartment.
- Black laminate base has access holes on both ends and back for cord and plug.
- · Natural Oak laminate only.
- Choice of (10) Pro-Form[™] vinyl colors.

MODEL	LENGTH	DEPTH	HEIGHT
4079	78"	30"	31"



B. MODEL A9064



C. MODEL A9065



COMMON FEATURES

SPLIT LEG TABLES – MODELS A9064, A9065, A9068:

 Extremely versatile and functional split leg design.



 Each padded leg rest can be adjusted and locked into desired angle up to 45° with spring-loaded lock-in lever.

600

Comfortable 2" thick high-density foam upholstered tops.

COLORS: (SEE INSIDE FRONT COVER)

Laminate - Natural Oak, Wild Cherry or Folkstone Gray

Upholstery – 10 Pro-Form[™] vinyl colors

OPTIONS

MODELS A9064, A9065, A9068 *69 - Gas Spring Backrest 0° to 75° (400 *99-L - Left hand table (A9068 only) *A946 - (1) Lock/Bolt for Double Doors

*FACTORY INSTALLED

A. MODEL A9068 Split Leg table

SPECIAL FEATURES:

Six hinged doors with interior shelf offer concealed supply storage.
Solid finished back.

MODEL	LENGTH	WIDTH	HEIGHT
A9068	78"	30"	31"

B. MODEL A9064 Split Leg Table

SPECIAL FEATURES:

• Six hinged doors with interior shelf on both sides of table.

MODEL	LENGTH	WIDTH	HEIGHT
A9064	78"	30"	31"

C. MODEL A9065 Split Leg Table

SPECIAL FEATURES:

• Open shelf storage, accessible from both sides of the table.

MODEL	LENGTH	WIDTH	HEIGHT
A9065	78"	30"	31"

SPLIT LEG TABLES

COMMON FEATURES

SPLIT LEG TABLES – MODELS A9063, A9088, A9718:

- Extremely versatile and functional split leg design.
- Each padded leg rest can be adjusted and locked into desired angle up to 45° with spring-loaded lock-in lever.
- Comfortable 2" thick high-density foam upholstered tops.

COLORS: (SEE INSIDE FRONT COVER)

Upholstery – 10 Pro-Form[™] vinyl colors

A. MODEL A9063 Split Leg Table



 Sliding doors on one side. Solid finished back.
 Natural Oak, Wild Cherry or Folkstone Gray laminate

iai i i i ate.			
MODEL	LENGTH	WIDTH	HEIGHT
A9063	78"	30"	31"

B. MODEL A9088 6-leg split leg table



- Six legs and full plywood subtop for added strength.
- Laminate shelf with (2) long stretchers.High-pressure laminate legs, apron and
- H-Brace.
- New enhanced leg hardware system.
- Legs have protective black vinyl scuff caps.
- Natural Oak or Wild Cherry laminate.

MODEL	LENGTH	WIDTH	HEIGHT
A9088	78"	30"	32"





PROTEAM

15

C. MODEL A9718

C. MODEL A9718 Electric Split Leg treatment table

- Durable motorized operation with "hands-free" foot control.
- Electro-hydraulic cylinders are built into (4) black powder coated steel legs with levelers.
- Fully enclosed oak laminate cabinet with sliding laminate doors.

MODEL	LENGTH	WIDTH	HEIGHT
A9718	78"	30"	27" - 37"

OPTIONS

MODELS A9063, A9088, A9718

- *69 Gas Spring Backrest 0° to 75°
- *99-L Left hand table (A9063 and A9718.)

*FACTORY INSTALLED



- · High-pressure laminate legs, apron, and H-Brace.
- · Laminate shelf with (2) long stretchers.
- New enhanced hardware system

•	Legs have	protective	black	vinyl	scuff	caps.	
---	-----------	------------	-------	-------	-------	-------	--

MODEL	LENGTH	WIDTH	HEIGHT
A9069	78″	30"	32"





B. MODEL A9070 KNEE AND HIP FLEXION SPLIT LEG OPEN STORAGE TABLE



- Open shelf storage, accessible from both sides
- of the table.

 Full plywood subtop for added strength. 					
MODEL	LENGTH	WIDTH	HEIGHT		
A9070	78"	30"	32"		



C. MODEL A9071 600 KNEE AND HIP FLEXION SPLIT LEG **CABINET TABLE**

· Six hinged doors with interior shelf offer concealed supply storage.

 Full subtop for increased strength. 					
MODEL	LENGTH	WIDTH	HEIGHT		
A9071	78″	30"	32"		

COMMON FEATURES

KNEE AND HIP FLEXION SPLIT LEG TABLES - MODELS A9069, A9070, A9071:

- · Extremely versatile and functional split leg design.
- Each padded leg rest can be adjusted and locked into desired angle up to 45° with spring-loaded lock-in lever.
- Each padded leg rest can be adjusted with hip and knee joints locked into 45° angle with easy to operate hinge mechanism.
- Full plywood subtop for added strength.
- 30" x 78" upholstered top with comfortable 2" thick high-density foam.

COLORS: (SEE INSIDE FRONT COVER)

Laminate - Natural Oak, Wild Cherry or Folkstone Gray (Folkstone - A9070, A9071 only) Upholstery – 10 Pro-Form[™] vinyl colors

3mm)

MODELS A9069, A9070, A9071 *69 – Gas Spring Backrest 0° to 75° / 🛺


*99-L - Left hand table (A9071 only) *A946 - (1) Lock/Bolt for Double Doors

*FACTORY INSTALLED

YES! HAUSMANN CAN DO THAT!

Hausmann's top priority is customer satisfaction and takes great pride in the ability to customize products to meet specific needs. Custom fabrication from a rough sketch, modifying existing products, sourcing high-pressure laminate or alternate color vinyl are routine requests handled expertly by Hausmann's team of engineers and product specialists. Hausmann can work with you, from start to finish, to make your concept for a product a reality.

Have an idea for a product? Yes, we can!




Design concept based on drawing submitted by the Packers. 3D models and renderings created by Hausmann.



Taping station being pre-assembled in the Hausmann factory.



Finished taping station with custom WilsonArt[®] Fonthill Pear laminate in the Packers athletic training room.





Football athletic training room shown with custom Nevamar® Washington Apple laminate and Boltaflex® Candy Apple vinyl.

- PROTEAM



Rose-Hulman Institute of Technology



Athletic training room shown with custom laminate side panel on end of taping station.

A. MODEL A9021 Recovery Table



- Foot end of table has a gas-spring elevating leg rest (up to 75°) and a full-extension modality drawer.
- Middle section of table has a side pullout drawer with (4) drying rods.
- Head section of table has (2) hinged doors with interior adjustable shelf offering concealed storage.
- · Finished laminate back.
- Table can be used as a traditional treatment table with elevating backrest by reversing the athlete.
- Grommet holes in foot section of table and in black laminate base allow for easy cord management.
- 2" thick high-density foam upholstered tops.
- 3mm high-impact PVC edges on all door and drawer fronts.

COLORS: (SEE INSIDE FRONT COVER) Laminate – Natural Oak and Wild Cherry

Upholstery – 10 Pro-Form™ vinyl colors

MODEL	LENGTH	WIDTH	HEIGHT
A9021	78″	30"	31"

OPTIONS

MODELS A9021, A9075, A9078

*A945 – Drawer Lock *A946 – (1) Lock/Bolt for Double Doors *99L – Left Hand Table

*FACTORY INSTALLED

B. PRO-OPTIONS TABLE

600 3mm

- Designed to allow you to customize to meet your specific needs... it is as easy as 1-2-3.
- Choose your cabinet style open cabinets on both sides of table or cabinets on one side of table with solid back.
- Add optional elevating split leg feature and/or gas-spring backrest.
- Customize each compartment with any of five possible options (some restrictions).
- 2" thick high-density foam uphostered tops.
- 3mm high-impact PVC edges on all doors and drawer fronts.
- COLORS: (SEE INSIDE FRONT COVER)

Laminate – Natural Oak and Wild Cherry Upholstery – 10 Pro-Form[™] vinvl colors

		5	
MODEL	LENGTH	WIDTH	HEIGHT
A9075	78″	30"	31"
A9078	78"	30"	31"

† NOTES (for A9075 Open Cabinet Table)

- Head section of table has offset support (not open pass-thru compartment – only option available is A962 Double Doors).
- If A965 or A966 Drawer option is selected, opposite side will be solid panel.

A. MODEL A9021

Scan QR code for

more information

B. MODEL A9075

Shown with options:

A962-(2) Doors, A963 - Pass thru drawer

STEP 1 - Choose your cabinet style


(standard configuration is flat top):

A9075 - Open Cabinet Both Sides[†]

• A9078 - Cabinet with Solid Back

A966 - (3) Drawers, 69 - Gas Spring Backrest, A961 - Split Leg

Shown as right hand table (standard configuration) Shown with option: 32 - Pillow



PROTE

G9 -Gas Spring Backrest, O° to 75°



ROTEA

A961 -Elevating Split Legs

STEP 3 – Customize your compartments with any of the following options (standard compartment configuration is open with (1) adjustable shelf)



- A962 Double Doors (includes (1) adjustable shelf)
- A963 Pass-Thru Modality Drawer (A9075 only – includes 2 grommets for electrical cord management)
- A964 Modality Drawer
 (A9078 only includes 2 grommets for
 electrical cord management)
- · A965 2-Drawers[†]
- A966 3 Drawers[†]

SPECIALTY TREATMENT TABLES



lables.			
MODEL	LENGTH	WIDTH	HEIGHT
4082	72″	24"	32"



- Fold-away design saves floor space.
- 2" high-density foam upholstered top.
- Sturdy, floor-supported wall mounting plate for added strength.
- (2) folding high pressure laminate legs.
- Projects only 8" when folded.
- Wall fasteners not included.
- Natural Oak laminate
- 10 Pro-Form[™] vinyl upholstery colors.

MODEL	LENGTH	WIDTH	HEIGHT
4075	72"	24"	31"

OPTIONS

- FOR MODELS 4082, 4075 *00 – Special Table Height
- (for 4082 only)
- ***04** Table Safety Strap
- *07 Nose Cutout (for 4075 only)
- 27 Paper Dispenser & Cutter Combo
- 32 Pillow, full size
- **35** Pillow, small

*FACTORY INSTALLED



Model 4075 (folded)

TAPING TABLES

A. MODEL 4061 TWO-PERSON TAPING TABLE / CONVERTIBLE TREATMENT TABLE

- · Solid (1) piece construction (not modular).
- (2) Padded taping sections 26" wide x 34" deep x 36" high with (1) removable middle top pad.
- Just add middle pad to convert taping table into a treatment table.
- (3) Storage sections with oak laminate interiors and shelf.
- (1) Center drawer on nylon rollers.
- Doors have concealed, adjustable, self-closing Euro-hinges.
- · Comfortable 2" thick high-density foam tops.
- 3mm high impact PVC edges on all door and drawer fronts.

COLORS: (SEE INSIDE FRONT COVER)

Laminate - Natural Oak

Upholstery - 10 Pro-Form™ vinyl colorsMODELLENGTHWIDTH*HEIGHT

4061	74"	34"	36"
*width of cabi	net is 30"		

OPTIONS

MODELS 4061, 4062

*A945 - Drawer Lock

*A946 - (1) Lock/Bolt for Double Doors

*FACTORY INSTALLED

B.MODEL 4062 Convertible taping Bench/ Treatment table



(3mm

- Multipurpose table functions both as a taping bench with backrest and as a treatment table.
- To use as a taping bench, raise the gas-spring padded backrest up to 75° with recessed
- control handles on each side of the table.Just lower the backrest and the table quickly
- converts into a full-length treatment table.Convenient front full extension pullout for
- taping supplies.
- Doors have concealed, adjustable, 110° self-closing Euro-hinges.
- (2) 21"W x 24"D open storage sections with adjustable shelves allow access from either side of the table.
- 2" thick high-density foam tops in (3) sections
 32"L seating section, 30"L elevating middle section, and 15"L head section.
- 3mm high impact PVC edges on all door and drawer fronts.

COLORS: (SEE INSIDE FRONT COVER)

Laminate - Natural Oak

Upholstery – 10 Pro-Form™ vinyl colors

MODEL	LENGTH	WIDTH	HEIGHT
4062	78"	27"	34"



Shown as Treatment Table with middle pad in place



S&W TREATMENT TABLES

COMMON FEATURES

S&W TREATMENT TABLES: · Natural wood finish is standard

on all S&W Tables. • H-Brace Tables have triple-bolted, solid

450

- hardwood legs.
- · Cabinet tables are quality hardwood plywood and veneer construction.
- Extra-firm 2" high-density foam tops with rounded corners for comfort.

COLORS: (SEE INSIDE FRONT COVER)

Upholstery – 10 Pro-Form[™] vinyl colors

MODEL		LENGTH	WIDTH	HEIGHT
ATT2472	10,	72″	24"	31"
ATT2478	<i>10</i> -	78"	24"	31"
ATT2772	<i>10</i> -	72"	27"	31"
ATT2778	10-	78"	27"	31"
ATT3072	10 ,	72"	30"	31"
ATT3078	1 0 7	78"	30"	31"

OPTIONS

- AFW Finished Wood Shelf
- AU Upholstered Shelf
- 27 Paper Dispenser/Cutter Combo
- *69 Gas Spring Backrest 0° to 75° (N/A on 24"W Tables)
- *A946 (1) Lock/Bolt for Double Doors

*FACTORY INSTALLED

A. MODEL ATT314

SPECIAL FEATURES:

- 600
- · Four door cabinet with center divider and one adjustable shelf on each side.
- · Leg levelers for uneven floors. Satin finish wire nulls

 Satin inish wi 	re pulls.			
MODEL	LENGTH	WIDTH	HEIGHT	
ATT314	72″	30"	31"	



400

B. MODELS ATT385, ATT386



- Includes 30" long Gas Spring
- Backrest, 0° to 75°, as standard feature.

•	Includ	des f	full-\	width	ιAU	up	hol	stere	d sł	nelf.

MODEL	LENGTH	WIDTH	HEIGHT
ATT385	72"	30"	31"
ATT386 (not shown) 78"	30"	31"



AFW Wood Shelf

AU Upholstered Shelf



B. MODEL ATT385 Shown with standard 69 - Gas Spring Backrest

S&W WHIRLPOOL TABLES



A. MODEL AT810 Whirlpool table



- Natural finish hardwood table with Oak HPL marine grade plywood top and steps.
- Water-resistant urethane finish.
- \cdot Leg levelers for uneven floors.
- 47" wide x 32" deep available for 22", 24" and 26" whirlpools.
- Please specify manufacturer, model number and size and height of tank when ordering. Custom sizes available.
 We add approximately 1 ¼" to the width of the
- tank as a clearance factor.

MODEL	LENGTH	WIDTH	HEIGHT	
AT810	47"	32"	Variable	
0/A footprint including stops is 56"L x /2"W				

O/A footprint including steps is 56"L x 42"W

MODEL AT815 Whirlpool table with cabinet



300

300

Same as AT810 but includes 2 shelf cabinet on backside of table as shown on Model AT825.
Shelf dimensions are 35"W x 12"H x 9"D.

MODEL	LENGTH	WIDTH	HEIGHT	
AT815 (not shown)	47"	32"	Variable	

O/A footprint including steps is 56"L x 42"W


MODEL AT820 Whirlpool table featuring Starboard® St Top

- Same as AT810 but with StarBoard® ST top and steps.
- StarBoard® ST is a high-density polyethylene sheet which is waterproof and scratch-resistant and features a textured surface.
 Standard color is Sanshade.

MODEL	LENGTH	WIDTH	HEIGHT		
AT820 (not shown)	47"	32"	Variable		
O/A for a transition of the state of a EC'' $V = (2'')A/$					

O/A footprint including steps is 56"L x 42"W

B. MODEL AT825 Whirlpool table featuring Starboard® St top with cabinet



- StarBoard® ST is a high-density polyethylene sheet which is waterproof and scratch-resistant and features a textured surface.
- Standard color is Sanshade as shown.Includes 2 shelf cabinet on backside of table.

MODEL	LENGTH	WIDTH	HEIGHT	
AT825	47"	32"	Variable	
O/A footprint including steps is 56"L x 42"W				

A. MODEL AT810			·I
>			
	>		
		TIL	
T	0	TIM	



B. MODEL AT825



POWER TREATMENT TABLES



COMMON FEATURES

POWER TABLES – MODELS 4718, 4719:

- Heavy-duty, motorized hydraulic Hi-Lo system with 10" lift stroke.
- Electrically activated hydraulic cylinders are built into (4) corner telescoping legs.
- Legs are black powder-coated and equipped with levelers for uneven floors.
- "Hands-free" foot control operation.
- Vinyl upholstered tops with 2" high-density urethane foam.
- Hospital grade plug.
- 110V, 60Hz, 2.4 amps.
- Available option: 27 Paper Dispenser & Cutter.

COLORS: (INSIDE FRONT COVER) Upholstery - 10 Pro-Form[™] vinyl colors.

A. MODEL 4719 HI-LO POWER PLINTH TABLE WITH SHELF

SPECIAL FEATURES:

- Bariatric 600 lb. lift capacity.
- Motorized table elevates from 27" to 37" high.
- Oak laminate apron with black laminate storage shelf.

MODEL	LENGTH	WIDTH	HEIGHT
4719	78"	30"	27"-37"

400

B. MODEL 4718 HI-LO "EXTRA LONG" TRAINERS TABLE

SPECIAL FEATURES:

- Motorized table elevates from 27" to 37" high.
- Extra-long 84" x 30" wide upholstered top is ideal for taller athletes.
- Fully enclosed Oak laminate cabinet with (4) drawers and sliding door storage.
- Optional #69 30"L Gas Spring Backrest.

MODEL	LENGTH	WIDTH	HEIGHT
4718	84"	30"	27"-37"





- Electric motorized height adjustment from 18 ½" (wheelchair height) to 35" high.
- "Hands-free" foot control operation.
- Adjustable 26" long headrest has dual gas spring operation and single lever activation.
- Headrest raises to 83°. Lowers to 40° below flat.
- Base has 4 ½" clearance for most patient hoists.
- Steel frame with powder coated cream finish.
- Stocking upholstery vinyl color is Dove Gray.
- Retractable locking casters with lever activation.
- Includes paper dispenser and cutter.

• 110V, 60 Hz, 3.15 Amps.

MODEL	LENGTH	WIDTH	HEIGHT
4710	78"	28"	18 ½" - 35"



22

POWERMATIC TREATMENT TABLES

COMMON FEATURES

POWERMATIC[®] TREATMENT TABLES --MODELS 4746, 4747, 4750, 4756, 4757 :



- The Powermatic[®] electro-mechanical drive mechanism provides an infinite number of height positions.
- Patients can mount or dismount without a footstep. You can work faster and more efficiently without back strain or fatigue.
- · 2" High-density urethane foam vinyl top.
- "Hands free" foot control operation.
- Includes Paper Dispenser and Cutter.
- Hospital grade plug.
- · 110V, 60Hz, 3.5 Amps. U.L. Listed.
- New contoured black base plate with
 (4) levelers and satin aluminum trim.

COLORS: (INSIDE FRONT COVER)

Laminate - Folkstone Gray 4746, 4747 only: choice of any high pressure laminate solid color or woodgrain

Upholstery – 10 Pro-Form[™] vinyl colors.

OPTIONS

- FOR MODELS 4746, 4747, 4750, 4756, 4757:
- *04 Table Safety Strap
- *07 Nose Cutout (4750 only)
- *15 (1) Safety Rail (4746, 4747 only) (Order [2] #15 if both sides)
- 32 Pillow, full size
- 35 Pillow, small
- *99L Left Hand Table Model

*FACTORY INSTALLED



A. MODEL 4746 Powermatic® treatment table With Gas-spring backrest

SPECIAL FEATURES:

- Gas-spring pneumatic 30" backrest elevates to 75° with new lever control on both sides.
- Power height adjusts from 26" to 38".
- · (2) laminate drawers on nylon rollers.
- Hinged door storage compartment.
- · Shipped fully assembled.

MODEL	LENGTH	WIDTH	HEIGHT
4746	76"	27"	26"-38"



 Same as Model #4746 except power 30" backrest can be elevated automatically to 75° with foot control.

MODEL	LENGTH	WIDTH	HEIGHT
4747	76"	27"	26"-38"

C. MODEL 4756



B. MODEL 4750 Powermatic® Table With Flat Top



SPECIAL FEATURES:

 Power height adjustment from 19" wheelchair height to 31" high.

MODEL	LENGTH	WIDTH	HEIGHT
4750	76"	27"	19"-31"

C. MODEL 4756 Powermatic® Table with GAS-Spring Backrest

SPECIAL FEATURES:

- Power height adjustment from 19" wheelchair height to 31".
- Gas-spring pneumatic 30" backrest elevates to 75° with new lever control on both sides.

MODEL	LENGTH	WIDTH	HEIGHT
4756	76"	27"	19"-31"

MODEL 4757 (not shown) POWERMATIC® TABLE WITH POWER BACKREST

 Same as Model 4756 except Power 30" backrest elevates to 75° with foot control

MODEL	LENGTH	WIDTH	HEIGHT	
4757	76″	27"	19"-31"	



450

COMMON FEATURES

HI-LO SECTIONAL TABLES -MODELS 6071, 6074:

- "Hands-free" foot control operation.
- · Open base allows closer access to patient.
- · Powder coated heavy duty steel frame
- Quiet and stable, single column actuator adjusts from 24" to 33" H.
- · Single column actuator saves space, no top shift when elevating.
- Hospital grade power cord.
- · 100V-240V, 50/60Hz, 3.15 Amps

A. MODEL 6071 2-SECTION HI-LO TREATMENT TABLE WITH ARMRESTS AND ADJUSTABLE GLIDES

B. MODEL 6074

SPECIAL FEATURES:

- · 2-Section Power Table can be used for a wide variety of purposes.
- · Folding armrests for extra patient comfort.
- Includes paper dispenser (factory installed).
- 2-Section top is 75" L x 28" W with 2" thick high density foam.
- Head/Backrest is 25" L x 28" W & elevates up to 67° with dual release friction lock.
- · Adjustable glides allow for easy leveling.
- · Agate Blue vinyl upholstery.

MODEL	LENGTH	WIDTH	HEIGHT
6071	75"	28"	24"-33"

B. MODEL 6074 3-SECTION MOBILE HI-LO TREATMENT TABLE WITH LOCKING CASTERS

SPECIAL FEATURES:

- 3-Section Power Table can be used for a wide variety of purposes.
- (4) 3" diameter, swiveling and locking casters.
- · 3-Section top is 78" L x 28" W with 2" thick high density foam.
- · Contoured face/nose cutout in head section.
- Head section, 19" L x 28" W adjusts from flat to 90° below horizontal.
- Middle section, 21" L x 28" W is stationary.
- Leg/Backrest section, 38" L x 28" W, elevates up to 70° with dual release friction lock.
- Agate Blue vinvl upholsterv.

MODEL	LENGTH	WIDTH	HEIGHT
6074	78"	28"	24"-33"



C. MODEL 6053 3-SECTION HI-LO MULTI-POSITION TREATMENT TABLE

- · Multi-position table can be used for a wide variety of therapy and treatment functions
- 76"L x 27"W 3-section 1 ½" thick high density gray vinyl top features contoured face/nose cutout in head section
- Motorized height adjustment from 19" wheelchair height to 37"H.
- "Hands-free" foot control operation.
- Adjustable head section, 15" x 27", raises up to 45° and lowers down to 90° below horizontal.



- 22"L x 27"W Middle seat section is stationarv
- 39"L x 27"W Leg/Backrest section raises easily up to 73° with sturdy friction locking device and dual release handles.
- Top sections have welded steel support frames.
- Includes (4) 3" diameter swiveling casters with locking brakes.
- 110V, 60 Hz, 1.8 amps. U.L. Listed.

MODEL	LENGTH	WIDTH	HEIGHT
6053	76"	27"	19"-37"

EXAMINATION TABLES

COMMON FEATURES

MAXIMUM VALUE EXAM TABLES – SERIES 4420, 4421, 4423, 4424;



- Removable vinyl top pad attached with VELCRO® brand fasteners.
- · Contoured padded head section.
- Backrest is pneumatic air-spring with dual adjustment controls from 0° to 75°.
- Built-in front safety tread footstep.
- Sliding legrest extension with removable pad attached with VELCRO® brand fasteners.
- Heavy duty concealed stirrups adjust in length and are laterally adjustable in (4) positions. Vinyl coated heelrests.
- Paper dispenser and paper cutter.
- Hospital grade duplex outlet. Note: LD models do not have outlet.
- Black laminate base with leveling legs.

COLORS: (SEE INSIDE FRONT COVER)

Laminate – Folkstone Gray Upholstery – 10 Pro-Form[™] vinyl colors

OPTIONS

FOR SERIES 4420, 4421, 4423, 4424

*15 - (1) Safety Rail (Order [2] if both sides)
*091 - Pelvic Tilt (Controls on both sides)
*097 - Stainless Steel Pan & Cutout

*FACTORY INSTALLED

A. SERIES 4423, 4424

SPECIAL FEATURES:

 (2) Drawers in front and (3) drawers on side as shown.

MODEL	LENGTH	WIDTH	HEIGHT
4423	55" - 71"	26"	33"
4423-LD	Same but less Duplex outlet		
4424	Left Hand Table		
4424-LD	Same but less Duplex outlet		

B. SERIES 4420, 4421

SPECIAL FEATURES:

• (2) Drawers in front and (1) door on side with shelf as shown.

MODEL	LENGTH	WIDTH	HEIGHT
4420	55"- 71"	26"	33"
4420-LD	Same but less Duplex outlet		
4421	Left Hand Table		
4421-LD	Same but less Duplex outlet		

Full-comfort pad with VELCRO® brand fasteners detaches easily for cleaning or replacement.





Pneumatic Air

Spring Backrest

A. MODEL 4423

B. MODEL 4420 Shown with option: 097 - Stainless Steel Pan

Upholstered, removable top pad

Slide-Out leg rest with pad, 21" wide by 16" deep

> 4-position laterally adjustable stirrups, concealed when not in use

Stainless steel irrigation pan (optional)

Storage cabinet with shelf

> Hospital grade duplex outlet (not supplied on "LD" models)

Built-In Safety-Tread Foot Step, 21" wide x 13" deep x 7¾ " high



A. SERIES 1440 Dual Lift Powermatic® Mat Platforms



- Power height adjustment from 19" wheelchair access height to 28"H.
- Two direct drive, synchronized lift mechanisms for unsurpassed power.
- Black laminate base with satin aluminum trim and levelers.
- Tear resistant, seamless nylon-reinforced vinyl upholstered top with 2" high-density urethane foam padding for comfort and rounded corners for safety.
- "Hands-free" foot control.
- · Hospital grade plug.
- 110V, 60Hz, 2.8 Amps. U.L. Listed

B. SERIES 1441 Econo electric hi-lo mat platforms



- Power height adjustment from 20" wheelchair height to 29"H.
- Two direct drive, synchronized lift mechanisms for greater power. 600 lbs. patient lift capacity.
- Tear resistant, seamless nylon-reinforced vinyl upholstered top with 2" high-density urethane foam padding for comfort and rounded corners for safety.
- · Sturdy black laminate base with leveling legs.
- "Hands-free" foot control.
- Hospital grade plug.
- · 110V, 60Hz, 1.3 Amps.

MODEL	WIDTH	LENGTH	HEIGHT
1441-46*	4'	6'	20"-29"
1441-47*	4'	7′	20"-29"
1441-57*	5'	7′	20"-29"

You must specify upholstery color with order.

MODEL	WIDTH	LENGTH	HEIGHT
1440-46-113	4'	6'	19"- 28"
1440-46-119	4'	6'	19"- 28"
1440-47-113	4'	7'	19"- 28"
1440-47-119	4'	7'	19"- 28"
1440-57-113	5'	7'	19"- 28"
1440-57-119	5'	7'	19"- 28"
1440-68-113	6'	8'	19"- 28"
1440-68-119	6'	8'	19"- 28"
Vou must speci	ify upholstor	v color wit	h order

You must specify upholstery color with order.



B. MODEL 1441-47-119



COMMON FEATURES MAT PLATFORMS -



- Extremely stable lift system elevates from 20" wheelchair height to 30" high.
- · "Hands-free" foot control.

SERIES 1427, 1428, 1429:

- · Seamless nylon-reinforced vinyl is wear and tear resistant.
- · 2" High-density urethane foam vinyl top.
- Steel frame and 4 corner legs have durable black powder coat finish.
- · Legs have levelers for uneven floors.
- · 4-Leg "open" design is patient sling-lift accessible on all four sides.
- Hospital grade plug.

COLORS: (SEE PAGE 24) Upholstery - 113 Blue, 119 Gray, 121 Green, 122 Black You must specify upholstery color with order.





A. SERIES 1427



SPECIAL FEATURES:

- · Electrically activated hydraulic cylinders are built into each (4) corner steel leg.
- · Motorized lift system can securely lift 750 lbs. from 20" wheelchair height to 30" high.
- "Hands-free" foot control operation.

MODEL	WIDTH	LENGTH	HEIGHT
1427-46*	4'	6'	20"-30"
1427-47*	4'	7'	20"-30"
1427-57*	5'	7'	20"-30"
1427-68*	6'	8'	20"-30"

B. SERIES 1428 ELECTRIC HI-LO MAT PLATFORM WITH POWERBACK

SPECIAL FEATURES:

• Electrically activated hydraulic cylinders built into each of 4 corner legs with 750 lbs. patient weight capacity.

750

• Wheelchair Accessible - Elevates smoothly from 21" wheelchair height to 31" H.

A. MODEL 1427-57-121

- Electric backrest can lift and elevate the back of an immobile patient. Backrest lift capacity is 400 lbs. Backrest adjusts from 1° to 70°.
- "Hands free" foot control operation for backrest and height adjustment.
- Ready for Use Ships fully assembled ready for use. Just Plug it in and Good to Go!

MODEL	WIDTH	LENGTH	HEIGHT
1428-47	4'	7′	21"-31"
1428-57	5'	7'	21"-31"



C. SERIES 1429 BARIATRIC 1100 LB. CAPACITY HEAVY-DUTY MAT PLATFORM



- · Electrically activated hydraulic cylinders are built into each (4) corner steel leg.
- · Motorized lift system can securely lift 1100 lbs. from 20" wheelchair height to 30" high.
- · Static load tested to 4400 lbs.
- · Seamless nylon reinforced vinyl is tear resistant.

 Plywood reinforced subtop for added strength. 				
MODEL	WIDTH	LENGTH	HEIGHT	
1429-57	5'	7'	20"-30"	
1429-68	6'	8'	20"-30"	

1100

COMMON FEATURES UPHOLSTERED MAT PLATFORMS -MODELS 1412, 1420, 1421:



- · 2" Thick high-density urethane foam.
- Upholstered in seamless heavy duty, . tear-resistant nylon-reinforced vinyl.
- High pressure laminate finish on all exterior surfaces.
- 18" height is convenient for therapist and wheelchair accessible for patient.

COLORS:

Laminate - 346 Natural Oak or 927 Folkstone Gray

Upholstery - 4 nylon-reinforced vinyl colors You must specify laminate and upholstery color with order.







A. SERIES 1412 UPHOLSTERED MAT PLATFORMS

SPECIAL FEATURES:

· All models have (6) legs for increased strength and stability (1412-46 has 4 legs).


MODEL	WIDTH	LENGTH	HEIGHT
1412-46	4'	6'	18"
1412-47	4'	7′	18"
1412-57	5'	7'	18"
1412-68	6'	8'	18"

B. MODELS 1420, 1421 SPACE SAVER MAT PLATFORM

SPECIAL FEATURES:

- · Fold-away design saves floor space.
- · Sturdy, floor-supported mounting plate.
- Three folding high pressure laminate legs.
- · Wall fasteners not included.

MODEL	WIDTH	LENGTH	HEIGHT
1420	4'	7'	18"
1421	5'	7'	18"
(projects 8" w	hen closed)		





C. MODELS 1418, 1419 **BARIATRIC SPACE-SAVER THERAPY PLATFORM**



- · Six high pressure laminate legs.
- · 2" thick high-density foam top.
- · 1000 lbs. patient weight capacity.
- · Laminate: 346 Natural Oak (shown) or 927 Folkstone Gray.
- Upholstery: Choice of (10) Pro-Form[™] Vinyl Colors.

MODEL	LENGTH	WIDTH	HEIGHT	
1418	78"	30"	18"	
1419	78"	36"	18"	



A. OPTION A931 PROPRINT™ GRAPHIC Shown on 1441-47

A. OPTIONS A931, A932, A933 Proprint™ graphic on mat platform

- Your team logo or graphics are digitally printed on mat vinyl using PROPRINT™ technology.
- Ultra-durable imaging uses premium ink and an unlimited color selection and gradient shading.
- Logo can be printed vertically or horizontally up to the size of the mat platform.
- Seamless nylon-reinforced vinyl is wear and tear resistant and available in 14 colors for logos on 4'x6' and 4'x7' mat platforms.
- Logos on 5'x7' and 6'x8' mat platforms are fully printed on white seamless nylon-reinforced vinyl; smaller mat platforms can also be fully printed using this same technology.
- Please note that fully printed mat platforms may be more susceptible to scratches and blemishes, particularly on the edges and at the corners.

OPTION DESCRIPTION

A931	Logo for 4'x6" or 4'x7' mat platform
A932	Logo for 5'x7' mat platform
A933	Logo for 6'x8' mat platform

PROPRINT™ NYLON REINFORCED VINYL COLOR CHOICES:

PP01	PP02	PP03	PP04	PP05	PP06	PP07
PP08	PP09	PP10	PP11	PP12	PP13	PP14

B. SERIES 1425 Crank Hydraulic Mat Platforms



- Durable, super-strength crank hydraulic system can lift 750 lbs. from 20" wheelchair height to 30" high.
- Easy-to-use hand crank height adjustment with fold away handle.
- Hydraulic cylinders are built into each (4) corner steel leg.
- · Seamless reinforced vinyl is tear resistant.
- Steel frame is durable black powder coated finish.
- · 2" High-density urethane foam top.
- · Legs have levelers for uneven floors.
- 4-Leg "open" design is patient sling-lift accessible on all four sides.

MODEL	WIDTH	LENGTH	HEIGHT
1425-46*	4'	6'	20"-30"
1425-47*	4'	7′	20"-30"
1425-57*	5'	7'	20"-30"
1425-68*	6'	8'	20"-30"

You must specify upholstery color with order.



SERIES 5580 Hugger® Gold-Line Wrist and Ankle Weights

- Vinyl weights color-coded and weight imprinted for ease of use.
- Extra long 17" strap with hook and loop fasteners allows you to double-strap weight in place for added security.
- Heavy-duty stitching assures greater strength and longevity.
- Extra-large brass hanging grommet.
- All weights 6 lbs. and over have (2) grommets
 & (2) straps with hook and loop fasteners and
 (6) vertical pockets.
- All weights are filled with smooth nonabsorptive steel pellets for comfort.
- Hugger® weights are hygienically sealed in individual plastic protective packages.
- Hugger® weights are sold individually or in sets.

MODEL 5580-75



Each weight is factory sealed. All weights are sold individually or in sets.

		ਿ		-
Je.	-			
19-	MODEL 5	580-1	_	

MODEL	LBS.	COLOR
5580-05	1/2	Tan
5580-1	1	Blue
5580-15	1½	Olive
5580-2	2	White
5580-25	2 1/2	Red
5580-3	3	Gold
5580-4	4	Turquoise
5580-5	5	Black
5580-6	6	Tan
5580-7	7	Blue
5580-75	7 ½	Orange
5580-8	8	Tan
5580-9	9	Green
5580-10	10	Brown
5580-12	12	Black
5580-150	15	Black
5580-200	20	Tan



MODEL 5581 (not shown) SET OF 16 HUGGER® WEIGHTS



Consists	of: (2) each of 1 lb., 1 ½ lb., 2 lb., 2 ½ lb.,
3 lb., 4 lb.	, 5 lb. and (1) each of 7 ½ lb. and 10 lb.
MODEL	ТҮРЕ
5581	Set of (16) Hugger® Weights

MODEL 5588 (not shown) SET OF 5 HUGGER® WEIGHTS

Consists of: (1) each of 1 lb., 2 lb., 3 lb., 4 lb., and

5 Ib.	
MODEL	ТҮРЕ
5588	Set of (5) Hugger® Weights

VINYL COATED DUMBBELLS

MODEL	LBS.	COLOR
5510-V	1	Pink
5511-V	2	Aqua
5512-V	3	Plum
5513-V	4	Grape
5514-V	5	Teal-Green
5515-V	6	Burnt Orange
5516-V	7	Bluish Green
5517-V	8	Burgundy
5518-V	9	Blue
5519-V	10	Ebony Black



Weight is embossed on each end

MODEL 5502 (not shown) SET OF 6 VINYL DUMBBELLS



MODELCONSISTS OF5502(2) each 1 lb. through 3 lb.



(1) each 1 lb. through 5 lb.

5503




MODEL 5505 SET OF 10 VINYL DUMBBELLS MODEL CONSISTS OF


5505 (1) each 1 lb. through 10 lb.



30

WEIGHTS, RACKS & DUMBBELL WAGON





A. MODEL 5555 ECONOMY WALL DUMBBELL RACK

10 lb. vinyl coated dumbbells.

easily to any wall.

B. MODEL 5555-100

CONSISTS OF

MODEL

5555

· Black powder coated metal wire rack attaches

(10) hooks are graduated in size to hold 1 to

53 ¼"

ACCESSORIZED WALL DUMBBELL RACK

(1) 5555 rack, (1) 5505 set of (10) Dumbbells

HEIGHT WIDTH

10"



DEPTH

3 1/2"

C. MODEL 5576 DUMBBELL WAGON



- · Oak laminate finish on all surfaces. • Upper 34 ½" x 22 ½" tray work surface
- accommodates up to (20) various size dumbbells.
- · Lower storage shelf with vertical divider.
- · Can store weights and accessories.
- · 3" diameter heavy duty casters.


MODEL	LENGTH	WIDTH	HEIGHT
5576	39"	24"	32"

MODEL 5576-100 ACCESSORIZED DUMBBELL WAGON

CONSISTS OF

E. MODEL 5545

5576 with (2) 5505 sets of (10) Dumbbells each



MODEL 5539 (not shown) **SET OF 3 EXERCISE WEIGHT BARS**



- · A series of 1 to 3 lb. easy-to-grip color-coded 36" long weighted bars.
- · Ideal for geriatric and nursing home use.
- Features built-in hanger end cap on one end . and anti-roll end cap on other end.

MODEL	CONSISTS OF	
5539	(1) each 1 lb., 2 lbs., 3 lbs.	

D. MODEL 5541

MODEL

5541

SET OF 5 EXERCISE WEIGHT BARS

- · A series of 1 to 5 lb. easy-to-grip color-coded 36" long weighted bars.
- · Ideal for rehabilitation exercise activity.

CONSISTS OF

• Features built-in hanger end cap on one end and anti-roll end cap on other end.

(1) each 1 lb., 2 lbs., 3 lbs., 4 lbs, 5 lbs.

= -1 -
123





- · Each 5" ball can be gripped with one hand.
- PVC exterior, anti-slip, color-coded, inflatable

with s	sealed	weigh	t pel	lets.

MODEL	CONSISTS OF
5545	(1) each 1.1 lbs., 2.2 lbs., 3.3 lbs.,
	4.4 lbs. 5.5 lbs., and 6.6 lbs.

WEIGHT RACKS / REP BANDS




A. MODEL 5570 **WEIGHT RACK**



- · Almond coated pegboard on both front and back side.
- · Natural Oak laminate base and frame.
- Includes (11) pegboard safety hooks. 2" diameter swiveling easy-rolling casters

MODEL WIDTH DEPTH HEIGH	
	Т





CONSISTS OF

(1) 5570 rack, (1) 5581 Set of (16) Hugger® Weights

COMMON FEATURES

LATEX-FREE RESISTIVE EXERCISE BANDS:



- · For multi-functional exercise use.
- Bands are 50 yds. long x 4" wide. Supplied . in dispensing box.
- Available in (5) color-coded resistance • levels.
- · Non-slip, non-latex elastomer.
- No powder, no odor.
- Use this clinically tested allergen-free safer alternative.





All weights and dumbbells shown must be ordered as accessories

B. MODEL 5565 "COMBO" WEIGHT/ **DUMBBELL MOBILE RACK**

hooks for (20) dumbbells.

Oak laminate frame & base.

(11) hooks.

MODEL

5565



MODEL 5565-100 ACCESSORIZED COMBO RACK MODEL CONSISTS OF



5565 Rack, (1) 5581 Set of (16) Hugger® Weights, (2) 5505 Sets of (10) Dumbbells each

MODEL 5571 PEGBOARD HOOKS

Set of (11) additional Pegboard Hooks with ball . tip ends for safety.

MODEL	LENGTH	QUANTITY
5571	4"	11

C. MODEL 5531

SET OF (5) REP BANDS™

Consists of the following:		
LEVEL	COLOR	
Extra-Light	Peach	
Light	Orange	
Medium	Green	
Heavy	Blue	
Extra-Heavy	Plum	

storage of optional Hugger® Weights, includes

Black laminate panel on back equipped with

WIDTH DEPTH HEIGHT

18"


Resistive Exercise Products

3" diameter swiveling, easy-rolling casters.

24"



57"



COMBINATION WEIGHT RACKS

COMMON FEATURES

COMBINATION WEIGHT RACKS MODELS 5558, 5560:

· Almond coated pegboard on front of rack for (16) optional Hugger® weights.

300

- (11) Pegboard safety hooks.
- (10) Dumbbell hooks on left side and (10) hooks on right side (total of 20 hooks). Dumbbells are optional.
- · Durable Black laminate finish.
- Base plate has bumper "T" molding on all four sides.
- · 3" diameter swiveling, easy-rolling casters.

A. MODEL 5560

MULTI-PURPOSE COMBINATION RACK

SPECIAL FEATURES:

- Built-in dispensers for #5531 set of (5) REP Bands™ (see page 32). Includes scissors.
- · Full view distortion-free glass mirror with ANSI "shatter stop" safety backing.

MODEL	WIDTH	DEPTH	HEIGHT
5560	34"	19"	71″

MODEL 5560-100 ACCESSORIZED MULTI-PURPOSE COMBINATION RACK

CONSISTS OF

(1) 5560 Rack, (1) 5581 Set of (16) Hugger® Weights, (1) 5531 Set of (5) REP Bands™, (2) 5505 Sets of (10) dumbbells each

B. MODEL 5558 COMBINATION RACK WITH WHITEBOARD

SPECIAL FEATURES:


- · HPL Whiteboard is a nonmagnetic surface and resists scratching.
- 48" high whiteboard offers a large blank space that's suitable for notes, diagrams, or therapy instructions. Shelf to hold markers and eraser. Built-in dispensers for #5531 set of (5) REP
- Bands[™] (see page 32). Includes scissors.

MODEL	WIDTH	DEPTH	HEIGHT
5558	34"	19"	71"

MODEL 5558-100 ACCESSORIZED RACK WITH WHITEBOARD

CONSISTS OF

(1) 5558 Rack, (1) 5581 Set of (16) Hugger® Weights, (1) 5531 Set of (5) REP Bands™, (2) 5505 Sets of (10) dumbbells each



Front view

Reverse view

All weights, dumbbells and bands must be ordered as accessories.



4-SIDED WEIGHT RACKS



A. MODEL 5557 Rotating Weight Rack

- Rotates 360° for ease of use.
- Built-in dispensers for (5) REP Bands[™].
- Hooks for (20) dumbbells.
- Pegboard with (6) hooks for weights.
- Stationary base plate has T-molding.

MODEL	WIDTH	DEPTH	HEIGHT
5557	20"	20"	56"

MODEL 5557-100 ACCESSORIZED RACK

(1) 5557 Rack, (2) Sets of 5588 Hugger® Weights,
 (2) 5505 Sets of (10) Dumbbells each, (1) 5531 Set of (5) REP Bands™

B. MODEL 5566 Econo Weight Rack

250



Pegboard with (11) hooks for Hugger® Weights.

- (20) Dumbbell hooks on reverse side of rack.
- (20) Dumbbell hooks on reverse side
- Built-in dispensers for (5) REP Bands™.
- Built-in compartment for weight bars.Black laminate finish. Base plate has
- T-molding.

• (4) 3" diameter swiveling casters.

()	5		
MODEL	LENGTH	DEPTH	HEIGHT

5566	30"	18"	59½"

All weights and dumbbells shown must be ordered as accessories

MODEL 5566-100 ACCESSORIZED RACK CONSISTS OF

10

(1) 5566 Rack, (1) 5581 Set of (16) Hugger®
 Weights, (2) 5505 Sets of (10) each 1 to 10 lb.
 dumbbells, (1) Set of 5531 Set of (5) REP™ Bands,
 (1) 5539 Set of (3) Weight Bars (1 lb, 2 lb, 3 lb)



C. MODEL 5592 Compact Weight Rack



• Pegboard with (11) hooks for Hugger® Weights.

- Built-in dispensers for (5) REP BandsTM.
- Hooks to hold (20) dumbbells.
- · Storage compartment for weight bars.
- Laminate shelves with (3) plastic storage tubs.
- Black laminate finish. Base plate has T-molding.
- (4) 3" diameter swiveling casters.

MODEL	LENGTH	DEPTH	HEIGHT
5592	32"	24"	68"

MODEL 5592-100 ACCESSORIZED RACK CONSISTS OF



(1) 5592 Rack, (1) 5581 Set of (16) Hugger®
 Weights, (2) 5505 Sets of (10) Dumbbells each,
 (1) 5531 Set of (5) REP Bands™, (1) 5541 Set of (5)
 Weight Bars, (1) 5545 Set of (6) "Soft-Ball" Hand
 Weights











A. MODEL 5562 Front view

Front view

MODEL 5562 Reverse view

All weights and dumbbells shown must be ordered as accessories





HHIDD

MODEL 5563 Reverse view

MODEL 5562-100 Accessorized Rack



CONSISTS OF

(1) 5562 Mega-Rack, (1) 5581 Set of (16) Hugger® Weights, (1) 5531 Set of (5) REP Bands™, (2) 5505 Sets of (10) Dumbbells ea., (1) 5539 Set of (3) Weight Bars

B. MODEL 5563 "Grand Stand" storage/weight rack

FEATURES:

- Full height 15" wide Oak laminate door with lock and self-closing Euro-hinges.
- Laminate cupboard with (5) adjustable laminate shelves.
- 3" high wire basket at top of rack holds miscellaneous supply items.
- MODEL
 LENGTH
 WIDTH
 HEIGHT

 5563
 32"
 32"
 72"*
- *75" Including Wire Basket

MODEL 5563-100 Accessorized Rack

CONSISTS OF

(1) 5563 Rack, (1) 5581 Set of (16) Hugger® Weights, (1) 5531 Set of (5) REP Bands™, (2) 5505 Sets of (10) Dumbbells ea., (1) 5539 Set of (3) Weight Bars



Front view

All weights and dumbbells shown must be ordered as accessories







- · Durable laminate finish.
- Top laminate surface is 36" W x 18"D.
- · 3mm PVC edging on all exposed edges.
- (5) High impact polystyrene 18" W drawers.
- · 18"W Oak laminate door with lock.
- Grab-bar and durable 4" locking casters for mobile use.
- · Hooks for (10) dumbbells.
- · Almond coated pegboard with (10) 2" safety hooks for Hugger® Weights.

MODEL	LENGTH	DEPTH	HEIGHT
5556	44"	24"	42"



MODEL 5556-100 ACCESSORIZED RACK

CONSISTS OF

(1) 5556 Rack, (2) 5503 Sets of (5) Dumbbells, (1) 5541 Set of (5) Weight Bars, (2) 5588 Sets of (5) Hugger® Weights.

300 A. MODEL 5569 MULTI-PURPOSE WEIGHT/ STORAGE RACK

- · Almond coated pegboard on front of rack for (16) optional Hugger® weights.
- (11) Pegboard safety hooks.
- (10) Dumbbell hooks on left side and (10) hooks on right side (total of 20 hooks). Dumbbells are optional.
- Durable Black laminate finish. .
- Base plate has bumper "T" molding on all four . sides.
- . 3" diameter swiveling, easy-rolling casters.
- · Two-door storage cabinet with lock.
- Doors have 165° self-closing Euro-hinges. .
- (2) Adjustable shelves and (1) fixed shelf.

MODEL	WIDTH	DEPTH	HEIGHT
5569	30"	24"	61"

MODEL 5569-100 ACCESSORIZED WEIGHT/ STORAGE RACK CONSISTS OF

(1) 5569 Rack, (1) 5581 Set of (16) Hugger® Weights, (2) 5505 Sets of (10) dumbbells each

Reverse view

(REP Resistive Exercise Bands™ must be ordered as optional accessory. See page 32.)



C. MODEL 6695 MOBILE CABINET FOR SPLINTING AND SUPPLIES

- · Grab-bar handle for mobility on left side.
- Door cabinet on right with lock.
- High impact polystyrene drawer. •
- (3) Vertical cubby sections for splinting material . behind hinged door.
- · Folkstone Gray laminate.
- · Durable dowel and glue construction.
- · 4" twin-wheel locking casters.
- Splinting pan is not included.

MODEL	LENGTH	DEPTH	HEIGHT
6695	42"	24"	31"

36



A. MODEL 8250

B. MODEL 8251

C. MODEL 8252

(3mm)

300

D. MODEL 8253

A. MODEL 8250 CABINET

SPECIAL FEATURES:

- Includes upper dispenser for optional REP Bands[™].
- (2) 10-hook metal wire racks on both interior sides to hold (20) optional dumbbells.
- Almond-coated pegboard panel on inside back to hold (16) optional Hugger® Weights. Includes (11) pegboard hooks.
- 3mm PVC high-impact edging on all door fronts.

MODEL	WIDTH	DEPTH	HEIGHT
8250	32"	22"	78"

MODEL 8250-100 Accessorized Cabinet

CONSISTS OF

(1) 8250, (2) 5505 Sets of (10) Dumbbells each
 (1 to 10 lbs), (1) 5581 Set of (16) Hugger[®] Weights,
 (1) 5531 Set of (5) REP Bands[™]

B. MODEL 8251 CABINET

SPECIAL FEATURES:

Includes (8) plastic storage tubs (12" x 14" x
 (" high) and adjustable laminate shelves

4 nigh) and adjustable faminate sherves.				
MODEL	WIDTH	DEPTH	HEIGHT	
8251	32"	22"	78″	

COMMON FEATURES

THERA-WALL™ THERAPY STORAGE SYSTEM

- MIX & MATCH Modular System. Order 1, 2, 3, 4 or more cabinets to fill any wall space.
- Each oak laminate cabinet module is 32" wide x 22" deep x 78" high and is equipped with (2) Oak laminate hinged locking doors.
- Modular cabinets must be secured to the wall and can be detached and moved if necessary. Wall fasteners not included.
- 3mm PVC high-impact edging on all door fronts.

C. MODEL 8252 CABINET

SPECIAL FEATURES:

- Includes storage bin for optional weight bars and space for canes and crutches.
- Almond pegboard back is equipped with (11) pegboard safety hooks for misc. items.

MODEL	WIDTH	DEPTH	HEIGHT
8252	32"	22"	78″

MODEL 8252-100 Accessorized Cabinet

CONSISTS OF

(1) 8252, (1) 5541 Set of (5) Weight Bars.


D. MODEL 8253 CABINET

SPECIAL FEATURES:


• (3) Laminate storage drawers behind doors.

 (3) Laminate adjustable shelves. 				
MODEL	WIDTH	DEPTH	HEIGHT	
8253	32"	22"	78″	

BOLSTERS & PILLOWS









C. POSITIONING PILLOWS

13

- Firm high-density urethane foam with vinyl sewn covering.
- Color for the positioning pillows is Medium Blue.

MO	DEL TYPE	LENGTH	WIDTH	HEIGHT
33	Cube	12"	8"	3"
34	Cube	12"	8"	6"
36	20° Wedge	18"	24"	6"
37	20° Wedge	32"	20"	12"
38	Cube w/40° edge	e 32"	10"	12"
39	Cube w/40° edge	e 32"	20"	12"

D. MODELS 41, 42, 43 "Anti-slip" Positioning Bolsters



- Firm high-density urethane foam with Medium Blue vinyl cover.
- Bolsters have "anti-slip" bottom surface that stays in place on mats or tables.

мо	DEL TYPE	LENGTH	WIDTH	HEIGHT
41	45° Wedge	16"	8"	8"
42	Half-round	17″	8"	7"
43	20° Wedge	32"	20"	12"

E. MODEL 3003 Multi-function twin wedge set



- Two "twin" high-density urethane foam wedges covered in Medium Blue vinyl.
- Each 20" x 20" x 10"H wedge has VELCRO® brand fasteners so wedges can be attached back-to-back or put together as a cube.
- Wedges can also be used separately for head, leg or knee support.

MODEL	WIDTH	LENGTH	HEIGHT
3003	20"	20"	10″



"CHOICE LINE™" AIR-LIFT STOOLS MODELS 2153, 2157:

- Hand operated air-spring pneumatic height adjustment.
- \cdot 16" diameter seat, 3 ½" thick foam.
- \cdot (5) 2" twin wheel black nylon casters.
- Durable 5-leg, black composite base.

COLORS: (SEE INSIDE FRONT COVER) Upholstery – 10 Pro-Form[™] vinyl colors

A. MODEL 2153 Economy Air-Lift Stool



250

250

SPECIAL FEATURES:

• Single lever height control handle.

MODEL	DIAMETER	HEIGHT	
2153	16"	20"- 27"	
2153-S	16"	17"- 22"	

B. MODEL 2157 Economy Air-Lift Stool



SPECIAL FEATURES:

• 360° circular ring control handle.

MODEL	DIAMETER	HEIGHT
2157	16"	20"- 27"
2157-S	16"	17"- 22"

C. MODEL 2116 Revolving stool

- $\cdot \;$ Chrome steel with round foot ring.
- $\cdot~$ (4) 2" twin wheel black nylon casters.
- Rotate stool top to adjust height.

· Choice of (10) Pro-Form™ vinyl colors.

MODEL	DIAMETER	HEIGHT	
2116	16"	18"- 24"	







information



250

6" Logo Stool Top

D. MODEL 2151 Heavy-duty Air-lift stool

D. MODEL 2151

FACTORY INSTALLED

INNS

-L

FOR MODELS 2116, 2151, 2153, 2157

Only available in PROFORM vinyl.

A900 - 6" plate for stool top debossing (required)

- Debossed Logo



- Heavy-duty pneumatic lift with "D" chrome control handle.
- 16" diameter seat, 3 ½" thick foam seat.
- Durable aluminum 5-leg base.
- (5) 2 ½" twin wheel gray nylon casters.
- Heavy-duty 400 lb. weight capacity.

Choice c	of (10) Pro-Form™ vin	yl colors.	
MODEL	DIAMETER	HEIGHT	
2151	16"	20"-24"	
			_



E. MODEL 2156-707



- Padded backrest, 10" high x 15" wide.
- Air-spring pneumatic height adjustment.
- 360° circular ring control handle.
- Durable 5-leg, black composite base.
- $\cdot~$ 16" diameter seat, 3 % " thick foam.
- (5) 2" twin wheel black nylon casters.

Black vinyl upholstery.

MODEL	DIAMETER	HEIGHT
2156-707	16"	20" - 27"

C. MODEL 2116





E. MODEL 2156-707



- Anti-slip treaded, extra-large 19" x 15" x 8"H platform with rubber angled leg tips.
- · MRI SAFE...non-ferrous and non-magnetic.

C. MODEL 2410 NESTED FOOTSTOOLS

MODEL	LENGTH	WIDTH	HEIGHT
2018	19"	15″	8"

· Consists of (4) steps which "nest" together for

· Durable Oak laminate wood construction with

20"

20" 18'

14"

12'

LENGTH WIDTH HEIGHT

16'

16"

14"

12"

9"

8"

8"

6"

4"

2"

400

500

· Nested footstools are useful for stepping, crawling or climbing activities.

easy space-saving storage.

non-slip matting.

MODEL

2410 (closed)

Stool Sizes:



B. MODEL 2038 BARIATRIC XL SAFE STEP WITH HANDRAIL

SPECIAL FEATURES:

- Weighs only 10 ½ lbs. for ease of use.
- Heavy duty tubular aluminum 15" x 19" x 8" high platform with unique anti-tip design.
- · 36" high sturdy steel safety handrail snaps in without tools.



E. MODEL 1310

B. MODEL 2038



D. MODEL 2010 BARIATRIC FOOTSTOOL · Durable 1" bent steel tubing in mainte-

- nance-free Silver Vein finish with safety tips. • 9" high non-skid "anti-slip" ribbed platform.
- · Legs are cross-braced to accommodate up to 500 lbs

500 165.			
MODEL	LENGTH	WIDTH	HEIGHT
2010	14"	10"	9"

E. MODEL 1310 BALANCE ACTIVITY PLATFORM

- · Adjustable height 3' stainless steel handrails.
- · Black high-pressure laminate with safety treads

MODEL	WIDTH	DEPTH	HEIGHT
1310	36"	36"	28"-41"

40

PROTEAM CARTS



COMMON FEATURES PROTEAM CARTS -SERIES A9014, A9015,

A9016, A9017;



- Maintenance-free laminated surfaces on all cabinet exteriors and interiors.
- Top laminate surface: 23" W x 18" D.
- · All cabinets have 3mm PVC edging
- High impact polystyrene drawer interiors with rounded inside corners.
- 4" satin finish wire drawer pulls with oak cabinet, or black nylon wire pulls with gray cabinet.
- 2 ½" twin wheel hooded swivel nylon casters.

COLORS: (SEE INSIDE FRONT COVER)

Laminate – 346 Natural Oak or 5904 Wild Cherry (A9014/A9015 Only). 927 Folkstone Gray & 5904 Wild Cherry available for contract/special orders

A. SERIES A9014 Proteam mega cart

SPECIAL FEATURES:

PROTEAM

• (1) 5 ½" drawer & (1) adjustable shelf.

MODEL	WIDTH	DEPTH	HEIGHT
A9014-346	24 ½"	18 ¾"	32 ½"
A9014-5904	24 ½"	18 ¾"	32 ½"

C. SERIES A9016 Proteam Monster Cart

SPECIAL FEATURES:

• (2) 5 1/2" drawers & (1) adjustable shelf.

MODEL	WIDTH	DEPTH	HEIGHT
A9016-346	24 ½"	18 ¾"	38"



B. SERIES A9015 PROTEAM MEGA CART *SPECIAL FEATURES:*



B. MODEL A9015-5904

(1) 5 ½" drawer, (2) hinged doors &
(1) adjustable shelf.

MODEL	WIDTH	DEPTH	HEIGHT
A9015-346	24 1⁄2"	18 ¾"	32 ½"
A9015-5904	24 1⁄2"	18 ¾"	32 ½"

D. SERIES A9017 Proteam Monster Cart

SPECIAL FEATURES:

•	(2) 5 ½" drawers, (2) hinged doors &
	(1) adjustable shelf.

MODEL	WIDTH	DEPTH	HEIGHT
A9017-346	24 1⁄2"	18 ¾"	38"



C. MODEL A9016-346



GAME READY™ CARTS / TREATMENT CABINET



A. MODEL A9022-346

COMMON FEATURES

GAME READY™ CARTS -A9022, A9023

- (3mm) 180 Maintenance-free laminated surfaces on all cabinet exteriors and interiors.
- All cabinets have 3mm high impact PVC edging for increased durability.
- Drawer features 1" opening for easy hose storage and high impact polystyrene interior with rounded inside corners for easy cleaning.
- · Grommets on top rear and back for easy electrical cord access.
- · 4" satin finish wire drawer pulls.
- 2 1/2" twin wheel hooded swivel nylon casters.

COLORS:

Laminate - 346 Natural Oak or 909 Black



A. MODEL A9022 GAME READY[™] CART

SPECIAL FEATURES:

- Designed to hold two Game Ready™ units.
- Top surface measures 17" width by 21 ½" depth.
- 4" drawer measures 15 ½" width by 14 ½" depth.
- . Bottom pull out features 3 drying racks for Game Ready™ sleeves.

MODEL	WIDTH	DEPTH	HEIGHT
A9022-346	18 1⁄2"	22 ¾"	38 ½"
A9022-909	18 ½"	22 ¾"	38 ½"



ANT ON FEADY

B. MODEL A9023 GAME READY™ MEGA CART



B. MODEL A9023-909

SPECIAL FEATURES:

- Designed to hold two Game Ready[™] units or one Game Ready™ unit and modality.
- Top surface measures 23" width by 21 ½" depth.
- 4" drawer measures 21 ½" width by 14 ½ depth.
- · Bottom pull out features 4 drying racks for Game Readv™ sleeves

Carrie Ready	5100705.		
MODEL	WIDTH	DEPTH	HEIGHT
A9023-346	24 ½"	22 ¾"	38 ½"
A9023-909	24 ½"	22 ¾"	38 ½"



C. MODEL 8728 TREATMENT CABINET					
• Sail White	nigh pressur	e laminat	e top.		
• (1) Laminat	e drawer on	nylon roll	ers.		
 All-laminat 	e interiors w	ith (1) adj	ustable shelf.		
Laminate - 346 Natural Oak or 927 Folkstone Gray.					
MODEL	WIDTH	DEPTH	HEIGHT		
8728-346	30"	16"	30"		
8728-927	30"	16"	30"		
				- I	

COMMON FEATURES

SMART CARTS™ & BEDSIDE CABINETS-MODELS 9018, 9029, 9038:



- Top laminate surface: 17" W x 17 1/2" D.
- · All cabinets have 3mm high impact edging. · Easy clean polystyrene drawer interiors with rounded inside corners.
- 4" satin finish wire drawer pulls with Oak cabinet, or Black nylon wire pulls with Gray cabinet.
- 2 ¹/₂" twin wheel hooded swivel nylon casters on all mobile models.
- · Ships fully assembled.

COLORS:

Laminate - Natural Oak or Folkstone Gray. Optional 5904 Wild Cherry & 7909 Fusion Maple (9038-20 Only)

10

3mm

40

D. MODEL 9210 MOBILE COMPUTER WORKSTATION

- · Portable and efficient mobile workstation with adjustable height top.
- Pneumatic Lift with 360° ring hand control.
- · 3mm edging for hi-impact & hi-durability.
- Durable 5-leg, black composite base with locking casters.
- · Non-rotating work surface.
- · 2" drawer with 3mm high impact edging.
- Plenty of space for portable electronics and . equipment.
- 8" height extender available Part 9211 . (height 33"-43")
- Laminate Natural Oak, Folkstone Gray, Wild Cherry

MODEL	WIDTH	DEPTH	HEIGHT
9210-346	22"	20"	25"-35"
9210-927	22"	20"	25"-35"
9210-5904	22"	20"	25"-35"



A. MODEL 9038-20-346



A. MODEL 9038-20 MOBILE CABINET

SPECIAL FEATURES:

(2) 5 ½" drawers & (1) 8 ½" deep drawer.

. ,	. ,		
MODEL	LENGTH	WIDTH	HEIGHT
9038-20-346	18 ½"	18 ½"	30"
9038-20-927	18 ½"	18 ½"	30"
9038-20-5904	18 ½"	18 ½"	30"
9038-20-7909 🌒	18 ½"	18 ½"	30"

B. MODEL 9029-24 **MOBILE CART**

SPECIAL FEATURES:

- · Accommodates EKG or compact modality equipment.
- · Available only in Folkstone Gray.

· Drop lid door & open storage.

		•	
MODEL	LENGTH	WIDTH	HEIGHT
9029-24-927	18 ½"	18 ½"	30"

B. MODEL 9029-24-927





Close up of 3mm edging.

CARTS & CABINETS NOW STANDARD WITH 3MM EDGING

- · High impact edging for increased product durability.
- · All edges are trimmed and buffed which eliminates sharp corners.



All finished edges are 3mm and will not peel.

MODEL 9018-20 **MOBILE CABINET** (not shown)

SPECIAL FEATURES:

• (1) 5 1/2" drawer & (1) hinged door.

Includes	(1) 11 ½" deep adjustable shelf.	

MODEL	LENGTH	WIDTH	HEIGHT
9018-20-346	18 ½"	18 ½"	30"
9018-20-927	18 ½"	18 ½"	30"

C. MODEL 9029-8 **MOBILE CABINET**

SPECIAL FEATURES:

- (1) 5 ½" drawer
- Open storage.

1 5			
MODEL	LENGTH	WIDTH	HEIGHT
9029-8-346	18 ½"	18 ½"	30"
9029-8-927	18 ½"	18 ½"	30"

C. MODEL 9029-8-927



INSTALLATIONS















STORAGE CABINETS / LOCKERS



3mm

COMMON FEATURES

TALL STORAGE CABINETS – MODELS 8154, 8248-84:

- All laminate dowel and glue wood construction. Shipped fully assembled.
- Easy to clean laminate on all exterior and interior surfaces.
- · Adjustable laminate shelves.
- · Satin finish wire pulls.
- Cabinets & lockers must be secured to the wall and can be detached and moved if necessary. Wall fasteners not included.

COLORS: (SEE INSIDE FRONT COVER) Laminate – 346 Natural Oak, 927 Folkstone Gray, 5904 Wild Cherry

C. MODELS A9590, A9594 Pro-lockers



- Dowel and glue oak laminate wood construction with 3mm high-impact edging.
- Foot locker with vented storage and hinged laminate seat (upholstered on A9594) with hasp for padlock.
- Double prong clothes hooks.
- · Chrome clothes rod on A9590 only.
- Cabinets & lockers must be secured to the wall and can be detached and moved if necessary. Wall fasteners not included.
- Can be manufactured to be ADA compliant.
 Quantity discount available for large orders

•	Quantity	discount availab	le lor large	orders.
		WIDTH		HEICHT

	WIDIII	DEFIN	HEIGHT
9590	32"	22"	78"
9594	18"	24"	78"
.9594	18"	24"	

OPTIONS

MODELS A9590, A9594 A9591 – Locking Security Door (for A9590) A9595 – Locking Security Door (for A9594) A. MODEL 8248-84-927

A. MODEL 8248-84 Double door storage cabinet

SPECIAL FEATURES:

- (2) doors and (8) adjustable shelves.
- (4) laminate drawers on nylon rollers.
- · Durable, 5 knuckle institutional hinges.
- Optional Lock (1 Door) or Lock-Bolt (2 Doors) available.

(, , ,				
MODEL	WIDTH	DEPTH	HEIGHT	
8248-84-346	48"	23"	76"	
8248-84-927	48"	23"	76"	
8248-84-5904	48"	23"	76″	



B. MODEL 8154-346

B. MODEL 8154 Single door storage cabinet

SPECIAL FEATURES:

- $\cdot~$ (1) door and (5) adjustable shelves.
- Durable, 5 knuckle institutional hinges.
- · Optional Lock available.

MODEL	WIDTH	DEPTH	HEIGHT
8154-346	24"	23"	76"
8154-927	24"	23"	76"
8154-5904	24"	23"	76"



C. MODELS A9590, A9594 Shown with optional locking

security door



CONTACT US ABOUT CUSTOM LOCKER CONFIGURATIONS

PRO-LINE™ MEDICAL CABINETS



COMMON FEATURES

PRO-LINE[™] PROFESSIONAL CABINETS:

- · A versatile line of modular all-laminate, dowel and glue cabinetry.
- · Neutral color all-laminate interiors.
- · All cupboards (except sink) have (1) adjustable laminate shelf.
- 4" Satin finish wire pulls.
- · Self-closing laminate drawer slides with nylon rollers and positive stops.

3mm

- Doors have concealed, self-closing 110° 3-Way adjustable Euro-Hinges.
- · All doors and drawer fronts have 3mm high impact edging.
- · All medical cabinets are "No Cancel/No Return."

COLORS: (SEE COLOR CHART THIS PAGE)


Cabinet - Folkstone Gray, Natural Oak or Wild Cherry laminate on all exposed surfaces including sides.

Countertop - Folkstone Gray square-edged top or choice of eight high pressure laminates with matching 3mm PVC edging.

920 ALMOND

966 BURGUNDY


D79 GREEN





IF SINGLE DOOR CABINET, YOU MUST SPECIFY DOOR HINGE "R" (RIGHT SHOWN) OR "L" (LEFT). ALL COUNTERTOPS HAVE 1/2" OVERHANG EACH END. SEND SKETCH AND DIMENSIONS WITH ORDER.

PRO-LINE™ MEDICAL CABINETS





CABINETRY INSTALLATIONS











Hausmann Industries is a premium certified member of the Architectural Woodwork Institute (AWI).

INDEX

MODEL #	PAGE	QS	NP/ND	MODEL #	PAGE	QS	NP/ND	MODEL #	PAGE	QS	NP/ND	MODEL #	PAGE	QS	NP/ND	MODEL #	PAGE	QS	NP/NI
30-43	38	3		5502	30	1		8250-100	37			A9023	42	10		A9531	4		
1310	40	3		5503	30	1		8251	37			A9024	10	5		A9532	4		
1412	28	5		5505	30	1		8252	37			A9041	10	5		A9533	4		
1418	28	10		5510-V to	70			8252-100	37			A9063	14	10		A9534	4		
1419	28	10	NP	5519-V	30	1		8253	37			A9064	13			A9535	4		
1420	28	5		5531	32	1		8728	42	3		A9065	13			A9540	4		
1421	28	5		5539	31	1		9018-20	43	5		A9066	11			A9541	4		
1425	29	10		5541	31	1		9029-8	43	5		A9067	11			A9542	4		
1427	27	10		5545	31	1		9029-24	43	5		A9068	13			A9543	4		
1428	27		NP	5555	31	1		9038-20	43	5		A9069	15		NP	A9544	4		
1429	27	10		5555-100	31	1		9210	43	10		A9070	15		NP	A9545	4		
1440	26	5		5556	36	10		A900	9			A9071	15		NP	A9590	45		
1441	26	10		5556-100	36	10		A901	9			A9073	11			A9594	45		
2010	40	1		5557	34	10		A902	9			A9075	17		NP	A9718	14		
2018	40	3		5557-100	34	10		A903	9			A9078	17		NP	AT810	21		
2038	40	3		5558	33	10	NP	A904-A908	9			A9087	10	10		AT815	21		
2116	39	3		5558-100	33	10	NP	A910	9			A9088	14	10		AT820	21		NP
2151	39	5		5560	33	5		A911	9			A9093	10	5	ND	AT825	21		NP
2153	39	3		5560-100	33	5		A912	9			A9097	10	5	ND	ATT2472	20	10	
2153-S	39	3		5562	35	5		A914	5			A9250	5			ATT2478	20	10	
2156-707	39	1		5562-100	35	5		A918	5			A9400	3		NP	ATT2772	20	10	
2157	39	3		5563	35			A923	5			A9401	3		NP	ATT2778	20	10	
2157-S	39	3		5563-100	35			A924	5			A9500	6			ATT3072	20	10	
2410	40	1		5565	32	1		A926	5			A9501	6			ATT3078	20	10	
3003	38	1		5565-100	32	1		A927	5			A9502	6			ATT314	20		
4061	19			5566	34	10		A929	5			A9503	6			ATT385	20	10	
4062	19		NP	5566-100	34	10		A931	29		NP	A9504	6			ATT386	20	10	
4075	18	5		5569	36	5		A932	29		NP	A9505	6			B-18 to BS-24	46		ND
4079	12		ND	5569-100	36	5		A933	29		NP	A9506	6			B-401	46		
4082	18	5		5570	32	1		A941	5			A9507	6			B-402	46		
4143	12	5		5570-100	32	1		A943	5			A9508	6			B-403	46		
4420	25	10		5571	32	1		A944	5			A9509	6			B-404	46		
4420-LD	25	10		5576	31	10		A945	5			A9510	6			GLR-A1	48		
4421	25	10		5576-100	31	10		A946	5			A9511	6			GLR-B4	48		
4421-LD	25	10		5580	30	1		A947	5			A9512	6			LOCK	47		
4423	25	10		5581	30	1		A961	17		NP	A9513	6			LOCK-BOLT	47		
4423-LD	25	10		5588	30	1		A962	17		NP	A9514	6			S	47		
4424	25	10		5592	34	10		A963	17		NP	A9515	6			S-LRG	47		
4424-LD	25	10		5592-100	34	10		A964	17		NP	A9516	6			SCS-24 to			
4710	22	3		6053	24	5		A965	17		NP	A9517	6			SCW-36	47		NP
4718	22			6071	24		NP	A966	17		NP	A9520	4			W-18 to W-24	47		
4719	22	10		6074	24	10	NP	A9002	10	5		A9520-346-24		5		LEGEND			
4746	23	-		6695	36	3		A9014	41	10		A9521	3,4					-10	YAQ
4747	23			7650-751	18	1		A9015	41	10		A9522	4			3 QUIC	K SHIP	- 3 [DAYS
4750	23	10		7650-752	18	1		A9016	41			A9523	4			_ •			
4756	23			8154	45	-		A9017	41			A9524	4			10 QUIC			
4757	23			8248-84	45			A9021	17		NP	A9525	4			NP NEW			
	20			52-10 04	1.5			, 13021	17		146	, (3323	т –						

©2016 HAUSMANN INDUSTRIES INC. The product dimensions and specifications in this catalog apply to the specific products and models described at the date of issue. Under our policy of continuous product improvement, all dimensions and specifications are subject to change without notice. Therefore, you should check on www.hausmann.com to ensure this catalog correctly describes the products currently available.



Tammy Boucher tammy.boucher@orthocanada.com 1.800.561.0310 x1700



CAT. #1016-PRO ©2016 HAUSMANN INDUSTRIES, INC. RECYCLED PAPER PRINTED IN THE USA



Scan for Savings













