



T635M TREADMILL

The T635M is rated for moderate cardiac and physical rehab and offers clients a low cost, value driven option when compared to our club quality product with many of the same functions and features.

- Bi-directional belt motion for retro walking and gait observations.
- Meets UL requirements for current leakage in medical settings
- Full length medical side rails with pediatric lower sections for better patient stability
- 400 lb user capacity, rated for 4-6 hours daily usage with 5 year parts, 3 year labor warranty
- Center handhold bars offers greater comfort for surgery patients

USER WEIGHT	DIMENSIONS	RUNNING SURFACE	INCLINE RANGE	SPEED RANGE	DRIVE MOTOR
400 lbs / 180 kg	83.5 x 39 x 57 in / 212 x 99 x 144 cm	22 x 61 in / 56 x 155 cm	0 to 15%	Reverse: 0.1 to 3.0 mph Forward: 0.1 to 10 mph	3.2 DC

Readouts

Calories, Speed, Time, Distance, Cal/HR, METs, Pace, Incline, Heart Rate, Weight Loss Zone, Cardio Zone, Active Zone

Workout Programs

Quick Start, Manual, Random (x 10,000), Interval (x3), Glute (x2) Hill (x3), 5K, 10K, Fat Burn, HRC (x3), Fitness Test (Bruce, Gerkin)

Features

- CSAFE port
- Medical Handrails
- Optional Wired Controller
- Forward and Reverse Speed Function (-3 to 10 mph)
- 0.1 SAFE-START™ Speed (forward & reverse)
- Contact and Telemetry HR
- 4-Ply 22 x 61 in. Belt
- USB Port (x2) for Device Charging
- Optional 3-Speed Fan
- ECO-GLIDE™ Auto Lubrication System
- Center Handle Bar with contact heart rate
- Optional SA Well +™
- Dual Headphone Jack with Built-in Volume Control

Warranty

Light Commercial ≤ 6 hrs/day
Frame Lifetime
MyFlex Plus™ Lifetime
Drive Motor 5 years

Parts 5 Years*
Wear Items 5 Years*
Labor 3 Years*
High Wear Items 90 days, no labor*

Tech Specs

Power Requirements 15A / 120v / 60Hz - Dedicated Circuit
Plug Requirements NEMA / 5-15P
Unit Weight 353 lbs / 160 kg
Step-Up Height 7.5 in / 19 cm
Micro Amp Leakage Components Medical Grade ≥ 300µA

Console



Tri-Color LED

Option



Wired Controller



*U.S. and Canada only.